

Kuruluş İzni Veriliş Tarihi: 30.12.2014		OXİGEN ANALİZ ÖZEL GIDA KONTROL LABORATUVARI FAALİYET KONULARI			Kuruluş İzni Yenileme Tarihi 03.08.2023	
Sicil No: 149						
Çalışma İzni Veriliş Tarihi: 20.01.2015		Adres: Çakmaklı Mah. Hadımköy Bağlantı Yolu Ufuk Plaza No:57 Kat:1 Daire:8 Büyüçekmece / İSTANBUL			Çalışma İzni Yenileme Tarihi: 06.05.2025	
Faaliyet İzin No:		Tel: (0212) 886 85 05-06 Fax: (0212) 886 85 04 E-Posta: info@oxigenanaliz.com Kep adresi: oxigenanaliz@hs01.kep.tr				
Sıra No	Analiz Adı	Metot Kodu	Metodun Kapsadığı Ürün / Ürün Grubu	Orijinal Metot Kaynağı / Kaynakları	Kullanılan Teknik *	Analiz Türü (Kalitatif/Kantitatif)
1	Nem (Rutubet) Tayini – Kuru Madde Tayini	FAL/SOP01/PR15	Hayvan Yemleri, Yemlik Yağlar, Yem Hammaddeleri ve Yem Katkı Madde Ürünleri	TS 6318, RG 29955 Yemlerin Resmi Kontrolü İçin Numune Alma ve Analiz Metotlarına Dair Yönetmelik	Gravimetrik	Kantitatif
2	GDO Tarama Analizi ¹ 2 (p35S, tNOS, pFMV)	MBAL/SOP01/PR15	Gıda ve Yem	Kit Metodu (Eurofins GeneScan GMO Screen RT 35S/NOS/FMV IPC Test Kit Manual), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
3	Bitki Spesifik Soya Geni Taraması ¹ 2	MBAL/SOP02/PR15	Gıda ve Yem	EURL Metot (QT-TAX- GM-002), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
4	Bitki Spesifik Mısır Geni Taraması ¹ 2	MBAL/SOP02/PR15	Gıda ve Yem	EURL Metot (QT-TAX- ZM-002), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
5	Bitki Spesifik Pamuk Geni Taraması ¹ 2	MBAL/SOP02/PR15	Gıda ve Yem	EURL Metot (QT-TAX- GH-015), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
6	Bitki Spesifik Şeker Pancarı Geni Taraması ¹ 2	MBAL/SOP02/PR15	Gıda ve Yem	EURL Metot (QT-TAX- BV-013), ISO 21569, ISO 24276	Real Time PCR	Kalitatif

7	Bitki Spesifik Kanola Geni Taraması ¹ 2	MBAL/SOP02/PR15	Gıda ve Yem	EURL Metot (QT-TAX-BN-012), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
8	Mezofilik Aerobik Bakteri Sayım Analizi	MAL/SOP01/PR15	Tüm Gıda, Gıda İle Temas Eden Her Türlü Yüzey ve Yem Maddeleri	TS EN ISO 4833-1 ISO 18593	Yayma / Dökme Plak	Kantitatif
9	Maya ve Küf Sayım Analizi	MAL/SOP02/PR15	Tüm Gıda ve Yem Maddeleri	TS ISO 21527-1 TS ISO 21527-2	Yayma / Dökme Plak	Kantitatif
10	<i>E.coli</i> Sayım Analizi	MAL/SOP03/PR15	Tüm Gıda ve Yem Maddeleri	ISO 16649-2	Yayma / Dökme Plak	Kantitatif
11	Koliform Analizi	MAL/SOP04-1/PR15	Tüm Gıda ve Yem Maddeleri	ISO 4831	EMS Yöntemi	Kantitatif
12	<i>E.coli</i> Sayım Analizi	MAL/SOP04-2/PR15	Tüm Gıda ve Yem Maddeleri	ISO 7251	EMS Yöntemi	Kantitatif
13	<i>E.coli</i> O157 Analizi	MAL/SOP05/PR15	Tüm Gıda ve Yem Maddeleri	ISO 16654	Yayma / Dökme Plak	Kalitatif
14	<i>Bacillus cereus</i> Analizi	MAL/SOP06/PR15	Tüm Gıda ve Yem Maddeleri	ISO 7932	Yayma / Dökme Plak	Kantitatif
15	<i>Clostridium perfringens</i> Analizi	MAL/SOP07/PR15	Tüm Gıda ve Yem Maddeleri	ISO 7937	Yayma / Dökme Plak	Kantitatif
16	<i>Cronobacter</i> spp. Aranması Analizi	MAL/SOP08/PR15	Süt ve Süt Ürünleri	ISO 22964	Yayma/Dökme Plak	Kalitatif
17	Enterobacteriaceae Familyasına Ait Mikroorganizmaların Sayım Analizi	MAL/SOP09/PR15	Tüm Gıda, Gıda İle Temas Eden Her Türlü Yüzey ve Yem Maddeleri	ISO 21528-2 ISO 18593	Yayma / Dökme Plak	Kantitatif
18	Koliform Bakteri Analizi	MAL/SOP10/PR15	Tüm Gıda, Gıda İle Temas Eden Her Türlü Yüzey ve Yem Maddeleri	ISO 4832 ISO 18593	Yayma / Dökme Plak	Kantitatif

19	Koagülaz Pozitif <i>Staphylococci</i> (<i>Staphylococcus aureus</i> ve diğer türleri) Analizi	MAL/SOP11/PR15	Tüm Gıda, Gıda İle Temas Eden Her Türlü Yüzey ve Yem Maddeleri	ISO 6888-1 ISO 18593	Yayma / Dökme Plak	Kantitatif
20	<i>Listeria spp</i> ve <i>Listeria</i> <i>monocytogenes</i> Analizi	MAL/SOP12/PR15	Tüm Gıda, Gıda İle Temas Eden Her Türlü Yüzey ve Yem Maddeleri	ISO 11290-1 ISO 18593	Yayma/Dökme Plak	Kalitatif
21	<i>Salmonella spp.</i> Analizi	MAL/SOP13/PR15	Tüm Gıda, Gıda İle Temas Eden Her Türlü Yüzey ve Yem Maddeleri	ISO 6579-1 ISO 18593	Yayma/Dökme Plak	Kalitatif
22	<i>Vibrio spp.</i> (<i>Vibrio</i> <i>cholerae</i> ve <i>Vibrio</i> <i>parahaemolyticus</i>) Analizi	MAL/SOP14/PR15	Tüm Gıda ve Yem Maddeleri	TS/ISO 21872-1	Yayma/Dökme Plak	Kalitatif
23	Termotolerant <i>Campylobacter</i> Türlerinin Analizi	MAL/SOP15/PR15	Tüm Gıda ve Yem Maddeleri	ISO 10272-1	Yayma/Dökme Plak	Kalitatif
24	Bitkisel Yağ Aranması	KAL/SOP02/PR15	Süt ve Süt Ürünleri	İşletme İçi Metot - “KAL/SOP02/PR15” (TS 7503'den modifiye edilmiştir.)	Kromatografik	Kantitatif
25	Sentetik Boya Miktarı Tayini [Sunset Yellow (E110), Azorubin, Tatrazine, Allura Red AC, Indigotine, Erythrosin, Brilliant Blue, Patent Blue V, Ponceau 4RC(E124), Kinolin, E- 123 Amarant]	KKMT/SOP06/PR15	Tüm Gıdalar	İşletme İçi Metot- ‘KKMT/SOP06//PR15 (NMKL 130, Waters Aplikasyon, Notu, Regnell'dan modifiye edilmiştir.)	Kromatografik	Kantitatif / Kalitatif
26	Boya Aranması	KKMT/SOP07/PR15	Tüm Gıdalar	Regnel 1976 Suda Eriyen Sentetik Renk Maddelerinin İzolasyon ve Teşhisi (Tarım Orman ve Köy İşleri Başkanlığı Koruma ve Kontrol Genel	Yün İpi	Kalitatif

				Müdürlüğü 1988.Gıda Maddeleri Muayene ve Analiz Metotları Bursa Sayfa 94-97)		
27	Biyojen Aminler (Histamin Tayini)	KKMT/SOP08/PR15	Su Ürünlerinde	Determination of histamine and biogenic amines in fish cubes implicated in a food borne poisoning Food Control 21 (2010)	Kromatografik	Kantitatif
28	Melamin Miktarı Tayini	KKMT/SOP09/PR15	Süt ve Süt Ürünleri, Yem ve Yem Katkı Maddesi, Bebek ve Küçük Çocuk Ek Gıdası	Journal Of AOAC international vol.91	Kromatografik	Kantitatif
29	Nitrat Tayini	KKMT/SOP10/PR15	Bebek ve Çocuk Ek Gıdaları, Yaş Sebze ve Meyve	NMKL165	Kromatografik	Kantitatif
30	Sorbik Asit ve Benzoik Asit Tayini	KKMT/SOP11/PR15	Tüm Gıdalar	NMKL 124.2	Kromatografik	Kantitatif
31	Kükürt Dioksit (SO ₂) Tayini	KKMT/SOP12/PR15	Tüm Gıdalar, Şarap ve Sirke, Bira (kuru soğan, pırasa, lahanası hariç)	AOAC 962.16	Titrimetrik	Kantitatif
32	Organoleptik Muayene	FAL/SOP02/PR15	Tüm Gıda ve Gıda Katkılarında, Tüm Yem ve Yem Katkılarında	TS EN ISO 5492	Organoleptik	Kantitatif
33	Laktik Asit Bakterilerinin Sayımı	MAL/SOP29/PR15	Gıda ve Yem	TS ISO 15214	Dökme Plak	Kantitatif
34	MON87701 Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP32/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-010) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
35	281-24-236 Pamuk Tip Belirleme Analizi ^{1 2}	MBAL/SOP34/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GH-001a) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
36	3006-210-23 Pamuk Tip Belirleme Analizi ^{1 2}	MBAL/SOP35/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GH-001b) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
37	Nitrit-Nitrat Analizi	KKMT/SOP16/PR15	Et ve Et Ürünler (Kürlenmiş Dil Hariç)	NMKL 165	Kromatografik	Kantitatif

38	Kafein Analizi	KKMT/SOP15/PR15	Gıdalar	İşletme içi metot KKMT/SOP15/PR15 (ISO 10727, Journal of AOAC International Vol.76, No:2)	Kromatografik	Kantitatif
39	Yapay Tatlandırıcı (Asesülfam-K, Sakarin, Aspartam) Tayini	KKMT/SOP18/PR15	Gıdalar	İşletme içi metot KKMT/SOP18/PR15 (TS EN 12856, Journal of AOAC International Vol.76, No:2)	Kromatografik	Kantitatif
40	Potasyum İyodat Analizi	KKMT/SOP17/PR15	Yemeklik Tuz	TS 933	Titrimetrik	Kantitatif
41	Natamisin Analizi	KKMT/SOP23/PR15	Süt Ürünleri, Kürlenmiş ve Kurutulmuş Et Ürünleri	İşletme içi metot KKMT/SOP23/PR15 (TS EN ISO 9233-2)	Kromatografik	Kantitatif
42	Asitlik Tayini	KAL/SOP03/PR15	Gıdalar	TS 2282 TS 1620 TS 522 TS 4500 TS 5000 TS 1880 EN 13188 TS 774 TS ISO 7305 TS 3036 TS 1330 TS 1329 TS 4680 TS 591 TS 1018 TS 1125 ISO 750 TS 7058 TS 1331 TS 1466	Titrimetrik	Kantitatif
43	Tuz Tayini	KAL/SOP05/PR15	Gıdalar	TS 774 TS 382 TS 3190 TS 2812 TS 5000	Titrimetrik	Kantitatif

				TS 591 TS 1466 TS EN ISO 5943 TS 2176 TS 1620 TS 1747-1 ISO 841-1 TS 11079 TS ISO 1738		
44	Tuz Tayini	KAL/SOP05/PR15	Yemler	21.01.2017 Tarih ve 29955 Sayılı Resmi Gazete 'de yayımlanan Yemlerin Resmî Kontrolü İçin Numune Alma ve Analiz Metotlarına Dair Yönetmelik	Titrimetrik	Kantitatif
45	Suda Çözünmeyen Madde Tayini	KAL/SOP07/PR15	Yemeklik Tuz	TS 1346 TS 933	Gravimetrik	Kantitatif
46	Yağ Tayini	KAL/SOP09/PR15	Kakao	TS 3076-1	Gravimetrik	Kantitatif
47	Yağ Tayini	KAL/SOP10/PR15	Gıda, Yem ve Yem Katkı Maddeleri	TS 7437 TS 1744 TS 6180 EN ISO 3947 TS 3190 TS 7800 RG 29955 Yemlerin Resmi Kontrolü İçin Numune Alma ve Analiz Metotlarına Dair Yönetmelik	Gravimetrik (Asit-Hidroliz Yöntemi)	Kantitatif
48	Yağ Tayini	KAL/SOP11/PR15	Gıda, Yem ve Yem Katkı Maddeleri	TS EN ISO 659 TS 1745 ISO 1444 TS 3713	Gravimetrik (Soxhelet)	Kantitatif
49	İyot Sayısı Tayini	KAL/SOP13/PR15	Hayvansal ve Bitkisel Yağlar	TS EN ISO 3961	Titrimetrik	Kantitatif
50	Sabun Miktarı Tayini	KAL/SOP15/PR15	Hayvansal, Bitkisel ve Yemlik Yağlar	TS 5038	Titrimetrik	Kantitatif

51	Sabunlaşma Sayısı Tayini	KAL/SOP16/PR15	Hayvansal, Bitkisel ve Yemlik Yağlar	TS EN ISO 3657	Titrimetrik	Kantitatif
52	Sabunlaşmayan Madde Tayini	KAL/SOP17/PR15	Hayvansal, Bitkisel ve Yemlik Yağlar	TS 4963	Gravimetrik	Kantitatif
53	Nişasta Tayini	KAL/SOP18/PR15	Et ve Et Ürünleri	TS 6812	Titrimetrik	Kantitatif
54	Ham Selüloz Tayini	KAL/SOP19/PR15	Gıda, Yem ve Yem Katkı Maddeleri	TS 4966, 21.01.2017 Tarih ve 29955 Sayılı Resmi Gazete’de yayımlanan Yemlerin Resmî Kontrolü İçin Numune Alma ve Analiz Metotlarına Dair Yönetmelik, TS 6932, TS ISO 15598, AOAC 962.09, TS EN ISO 6865, TS 324, TS 3076-1	Gravimetrik	Kantitatif
55	Protein Tayini	KAL/SOP21/PR15	Gıda, Yem ve Yem Katkı Maddeleri	TS 1620 TS 1748 ISO 937 TS EN ISO 8968-1	Titrimetrik (Kjheldal)	Kantitatif
56	Jelatin Aranması	KAL/SOP22/PR15	Süt ve Süt Ürünleri	AOAC 920.106	Organoleptik	Kalitatif
57	Toplam Uçucu Bazik Azot (TVB-N) Tayini	KAL/SOP23/PR15	Su Ürünleri	Balıkçılık Ürünlerine Ait Duyusal Özellikler ve Toplam Uçucu Bazik Azot Limitleri Tebliği (Tebliğ No: 2012/73)	Titrimetrik, Distilasyon Cihazı	Kantitatif
58	Sodyum Karbonat Aranması	KAL/SOP24/PR15	Süt	TS 1018	Organoleptik	Kalitatif
59	Dithiocarbamat Grubu Pestisit Analizi (Propineb)	KKMT/SOP22/PR15	Bebek Formülleri, Devam Formülleri ve Bebek ve Küçük Çocuk Ek Gıdaları	Validation of the method for the determination of dithiocarmates and rtiumam disulphide	Kromatografik	Kantitatif
60	Rutubet (Nem) (Uçucu Madde, Kuru Madde,	FAL/SOP03/PR15	Gıdalar	TS 7800 TS 1331 TS 3522	Gravimetrik	Kantitatif

	Katı Madde ve Kurutma Kaybı) Tayini			TS 1743 ISO 1442 TS 1018 TS 6800 TS 1561 TS 1562 TS 3076-1 TS ISO 7513 TS 1330 TS 7420 TS ISO 771 TS 1072 TS 1275 TS 1276 TS 1277 TS 1278 TS 1279 TS 1280 TS 3075 TS ISO 6479 TS EN ISO 5534 TS EN ISO 712 TS 5000 TS 1252 EN ISO 1666 TS 3190 TS 933 TS 2314 AOAC 979.12 TS 591 TS 12649 TS 3272 TS ISO 3728 TS 1329 TS 7058 TS 5389 TS EN ISO 665 TS ISO 3728		
61	Yoğunluk Analizi	FAL/SOP04/PR15	Süt Ürünleri, Alkollü İçkiler, Sıvı Yağlar, Meyve ve Sebze Ürünleri	TS 4959 TS 522 TS 5555	Gravimetrik	Kantitatif

				TS 342 TS EN 1131 TS 1018		
62	HCl de Çözünmeyen Kül Tayini	FAL/SOP08/PR15	Gıda ve Yemler	TS 2133 ISO 930 TS 13423 TS ISO 735 TS 2383 TS 5000 TS 1566 ISO 1577 TS ISO 5985 RG 29955 Yemlerin Resmi Kontrolü İçin Numune Alma ve Analiz Metotlarına Dair Yönetmelik TS 3076-1 TS ISO 763 TS 1466	Gravimetrik	Kantitatif
63	Kül Tayini	FAL/SOP05/PR15	Gıda ve Yemler	TS 2131 ISO 928 TS 1564 TS 1746 ISO 936 TS ISO 5984 RG 29955 Yemlerin Resmi Kontrolü İçin Numune Alma ve Analiz Metotlarına Dair Yönetmelik TS EN ISO 6498 TS ISO 6884 TS 13423 TS 3076-1 TS EN 1135 TS 522 TS 3792 TS ISO 749 TS 324 TS 5389 TS EN ISO 2171	Gravimetrik	Kantitatif

				TS 5816 EN ISO 3593 TS 1252 EN ISO 1666 NMKL 173, TS 2590, TS 2589		
64	pH Tayini	FAL/SOP06/PR15	Gıdalar	TS 2259 TS 5389 TS 3136 ISO 2917 TS 3792 TS 591 TS 2176 TS 1728 ISO 1842 TS 774 TS EN 1132	Potansiyometrik	Kalitatif
65	Refraktometrik Kırılma İndisi Tayini	FAL/SOP07/PR15	Hayvansal ve Bitkisel Katı-Sıvı Yağlar, Eteri Yağlar, Yemeklik Bitkisel Yağlar	TS EN ISO 6320, TS ISO 280 TS 894	Refraktometrik	Kantitatif
66	Uçucu Yağ Analizi	FAL/SOP09/PR15	Baharatlar, Çeşniler ve Tıbbi Bitkiler	TS EN ISO 6571	Gravimetrik	Kantitatif
67	Yabancı Madde Miktarı Tayini	FAL/SOP11/PR15	Tüm Gıdalar ve Yemler	TS 2066 TS 1279 TS 3689 ISO 7702 TS 4078 TS 12714 TS 2947 EN ISO 658 TS 1278 TS 1277 TS 5230 TGK Baharat Tebliği (Tebliğ No:2022/7) TS 3036 TGK Bebek Ve Küçük Çocuk Ek Gıdaları Tebliği TGK İnsan Tüketimine Sunulan Şekerlerin Analiz Metodları Tebliği TS 3882 TS 3706 TS 1620	Gravimetrik	Kantitatif

				TS 3415 TS 3267 TS 4600 ISO 3720 TS 1565 TS 3117 TS 3410 TGK Çeltik Tebliği TS EN ISO 927 TS 5389 TS 1466 TS 5000 TS 3687 ISO 7703 TS 1204 TS 3075 TS 3074 TS 312 TS 3223 TS 542 TS 10929 TS 5086 TS 1275 TS 2032 TS 13375 TS 3786 TS 314 TS 313 TS 4670 TS ISO 6658 TS 310 TS 309 TGK Kahve Ve Kahve Ekstraları Tebliği (Tebliğ No: 2016/7) TS 13423 TS 1315 TS ISO 2173 TS 2974 TS 2383 TS 3268		
--	--	--	--	---	--	--

				TGK Buğday Unu Tebliği (Tebliğ No: 2013/9) TS 10580 TS 141 TS 143 TS EN 12143 TS 3958 TS 3792 TS 10432 TS 2639 TS 3691 TS 1129 ISO 1026 TS 7896 TS 774 TS 1314 TS 2033 TS 142 TS 308 TS 316 TS 322 TS 487 TS 8596 TS 8597 TS 9278 TS 9696 TS 10991 TGK Pirinç Tebliği (2010/60) TS 2970 TS 933 TS EN ISO 3093 TS 2290 ISO 959-1 TGK Sofra Ve Gıda Sanayi Tuz Tebliği (Tebliğ No: 2013/48) TS 321 TS 9779 TS 522		
--	--	--	--	--	--	--

68	DP305423-1 Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP43/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-008) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
69	CV127 Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP44/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-011) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
70	MON87708 Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP46/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-012) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
71	MON87769 Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP47/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-002) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
72	GHB614 Pamuk Tip Belirleme Analizi ^{1 2}	MBAL/SOP48/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GH-006) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
73	Bitki Spesifik Pirinç Geni Taraması ^{1 2}	MBAL/SOP02/PR15	Gıda ve Yem	EURL Metot (QT-TAX-OS-017) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
74	GDO Tarama Analizi (p35S, tNOS,pFMV) ^{1 2}	MBAL/SOP56/PR15	Gıda ve Yem	Kit Metodu (Bosphore® GMO Screening Kit S35-TNOS- FMV), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
75	Enterobacteriaceae Sayımı	MAL/SOP34/PR15	Gıda ve Yem	ISO 21528-1	EMS Metodu	Kantitatif
76	<i>Pseudomonas</i> spp. Sayımı	MAL/SOP35/PR15	Et ve Et Ürünleri	TS EN ISO 13720	Yayma Plak	Kantitatif
77	Küf Sayımı	MAL/SOP36/PR15	Domates ve Domates Ürünleri (Domates suyu, domates püresi ve salça vb.)	AOAC 984.29	Mikroskopik (Howard Lamı ile)	Kantitatif
78	<i>Shigella</i> spp. Aranması	MAL/SOP43/PR15	Gıda ve Yem	ISO 21567	Yayma/Dökme Plak	Kalitatif
79	Rope Spor Sayımı	MAL/SOP44-2/PR15	Ekmek	TS 5000	EMS Yöntemi	Kantitatif

80	Rope Spor Sayımı	MAL/SOP44-1/PR20	Ekmek Mayası	TS 3522	EMS Yöntemi	Kantitatif
81	Karakteristik Mikroorganizma Sayımı	MAL/SOP46/PR15	Yoğurt	TS ISO 7889	Yayma / Dökme Plak	Kantitatif
82	Somatik Hücre Sayımı	MAL/SOP47/PR15	Çiğ Süt	TS ISO 13366-1	Mikroskopik	Kantitatif
83	Konservelerin İnkübasyon Testi (Ticari Sterilite Kontrolü)	MAL/SOP49/PR15	Konserve Ürünleri	TS 10524, Hayvansal Gıdalar İçin Özel Hijyen Kuralları Yönetmeliği	Yayma, Dökme, Tüp Metodu	Kalitatif
84	Enterococcus Sayımı	MAL/SOP55/PR15	Gıda ve Yem	NMKL 68	Yayma / Dökme Plak	Kantitatif
85	Maya ve Küf Sayımı	MAL/SOP56/PR15	Gıda, Çevre Yüzey Örnekleri (Gıda ile temas eden her türlü yüzey)	AOAC 2014.05, 3M Petrifilm Rapid Yeast and Mold Count Plate	Yayma Plak (Hızlı Metot)	Kantitatif
86	Tane/kg Sayısı - İrilik Derecesi Tayini	FAL/SOP12/PR15	Zeytin	TS 774, Türk Gıda Kodeksi Sofralık Zeytin Tebliği (Tebliğ No: 2014/33)	Gravimetrik	Kantitatif
87	1000 Tane Kütlesi Tayini	FAL/SOP13/PR15	Tahıl ve Baklagiller	TS EN ISO 520	Gravimetrik	Kantitatif
88	Boylama ve Elek Analizi	FAL/SOP14/PR15	Tüm Gıdalar	TS 933, TS 4500, TS 2970, Türk Gıda Kodeksi Pirinç Tebliği (Tebliğ No: 2010/60), TS 4201, TS 2640, TS 4078, TS 2974, TS 141, TS 142, Türk Gıda Kodeksi İrmik Tebliği(Tebliğ No: 2002/21), TS 4681, Türk Gıda Kodeksi Bulgur Tebliği (Tebliğ No: 2010/34), Türk Gıda Kodeksi Baharat Tebliği (Tebliğ No: 2013/12), TS 382, TS 3076-1	Gravimetrik	Kantitatif

89	Briks (Suda Çözünebilen Katı Madde) Tayini	FAL/SOP15/PR15	Tüm Gıdalar	TS ISO 2173, TS 1466, AOAC 932.14 TS EN 12143	Refraktometrik	Kantitatif
90	Et Oranı Tayini	FAL/SOP16/PR15	Zeytin	TS 12343, TS 12352	Gravimetrik	Kantitatif
91	Çemen Oranı Tayini	FAL/SOP17/PR15	Pastırma	TS 1071	Gravimetrik	Kantitatif
92	Elektrik İletkenliği Tayini	FAL/SOP18/PR15	Bal, Su, Şeker	TS 13366, TS 861, TS 266, TS 9748 EN 27888	Kondüktometrik	Kantitatif
93	Hektolitre Ağırlığı Tayini	FAL/SOP19/PR15	Tahıllar	TS EN ISO 7971-1, TS EN ISO 7971-2, TS EN ISO 7971-3	Gravimetrik	Kantitatif
94	Asitte Çözünmeyen Madde Miktarı Tayini	FAL/SOP20/PR15	Tuz	TS 1346	Gravimetrik	Kantitatif
95	Meyve ve Kaymak Oranı Tayini	FAL/SOP21/PR15	Şekerli Ürünler	TS 8444, TS 7800, TS 4186, TS 3958, TS 5230, TS 7780, TS 13555, Türk Gıda Kodeksi Kakao Ve Çikolata Ürünleri Tebliği (Tebliğ No: 2017/29), Türk Gıda Kodeksi Lokum Tebliği (Tebliğ No:2013/55)	Gravimetrik	Kantitatif
96	Kutu Doldurma Oranı Tayini	FAL/SOP22/PR15	Konserve Gıdalar	TS 2664	Gravimetrik	Kantitatif
97	Renk Tayini	FAL/SOP23/PR15	Bitkisel ve Hayvansal Yağlar, Bira	TS 2259, Lovibond PFX880 manual English v1.4	Spektrofotometrik	Kantitatif
98	Rutubet (Nem) Tayini	FAL/SOP24/PR15	Bal	TS 13365	Refraktometrik	Kantitatif
99	Suya Geçen Madde Miktarı Tayini	FAL/SOP25/PR15	Makarna, Bulgur	TS 1620	Gravimetrik	Kantitatif
100	Çözünbilme Oranı Tayini	FAL/SOP26/PR15	Süt Tozu, Çözünbilir Kahve	TS 1329, TS 5389	Gravimetrik	Kantitatif

101	Suda Çözünmeyen Madde Miktarı Tayini	FAL/SOP27/PR15	Tuz, Bal	TS 1346, TS 3036	Gravimetrik	Kantitatif
102	Süzme Ağırlığı Tayini	FAL/SOP28/PR15	Konserve Gıdalar	TS 382, TS 2664	Gravimetrik	Kantitatif
103	Su Ekstraktı Tayini	FAL/SOP29/PR15	Siyah Çay	TS ISO 9768	Gravimetrik	Kantitatif
104	Yağsız Katı Madde Miktarı Tayini	FAL/SOP32/PR15	Tereyağı	TS 1331	Gravimetrik	Kantitatif
105	Yaş ve Kuru Gluten Tayini	FAL/SOP33/PR15	Buğday Unu	TS EN ISO 21415- 2, TS EN ISO 21415- 4	Gravimetrik	Kantitatif
106	Erime (Kayma) Noktası Tayini	FAL/SOP34/PR15	Bitkisel ve Hayvansal Katı Yağlar	TS 2812, TS 11725, TS EN ISO 6321	Organoleptik	Kantitatif
107	Donmuş Etlerin Çözünmesiyle Ortaya Çıkan Su Miktarı Tayini	FAL/SOP35/PR15	Et	TS 8057	Gravimetrik	Kantitatif
108	Tortu Tayini	FAL/SOP36/PR15	Şarap, Yağlar	TS 894, TS 522	Organoleptik	Kalitatif
109	Külde Alkalilik Tayini	FAL/SOP37/PR15	Kakao Öğütülmüş	TS 3076-1	Gravimetrik	Kantitatif
110	Suda Çözünen ve Çözünmeyen Kül Tayini	FAL/SOP38/PR15	Çay, Kahve	TS 1565, TS 13423	Gravimetrik	Kantitatif
111	Suda Çözünen Külde Alkalilik Tayini	FAL/SOP39/PR15	Çay, Kahve	TS 1567, TS 13423, TS 5389	Gravimetrik	Kantitatif
112	Sülfatlandırılmış Kül Tayini	FAL/SOP40/PR15	Nişasta ve Nişastalı Ürünler, Gıda ve Gıda Katkı Maddeleri, Şeker ve Şekerli Ürünler, Şekerlemeler, Şurup, Laktoz ve Dekstroz	TGK İnsan Tüketimine Sunulan Şekerlerin Analiz Yöntemleri Tebliği (Tebliğ No:2002/26), TS 6361 EN ISO 5809	Gravimetrik	Kantitatif
113	Rutubet (Nem) Tayini – Kuru Madde Tayini	FAL/SOP41/PR15	Baharat ve Çeşni Veren Maddeler	TS ISO 939	Volumetrik	Kantitatif

114	Çeşni Maddelerinin Kütlesi-Çikolata Kısmı Tayini	FAL/SOP42/PR15	Bisküvi, Çikolata	TS 2383, TS 7800	Gravimetrik	Kantitatif
115	Sıvıda Su ve Yağ Oranı Tayini	FAL/SOP43/PR15	Balık Konserveleri	TS 353	Gravimetrik	Kantitatif
116	Hacim Genişlemesi Tayini	FAL/SOP44/PR15	Dondurma	TS 4265	Gravimetrik	Kantitatif
117	Fiziksel Kusurlar ve Özellikler Tayini	FAL/SOP45/PR15	Tüm Gıdalar	TS 4078, TS 12714, TS 309, TS 2974, TS 2284, TS 1276, TS 1275, TS ISO 2451, TS ISO 2451, TS 3411, 14.11 TS 3997, TS 1466, TS 11999, TS 3075, TS 2290 ISO 959-1, TS 3268, TS 141, TS 143, TS 3479 ISO 2591-1, TS 3415, TS 142, TS 3628, TS EN ISO 7301, TS 308, TS 3200, İşletme İçi Metot “FAL/SOP45/PR15” (TGK Tebliğleri, Ürün TS'leri, Ürün Spektleri)	Gravimetrik	Kantitatif
118	Ağırlık Kontrolü	FAL/SOP46/PR15	Yumurta	Türk Gıda Kodeksi Yumurta Tebliği (Tebliğ No: 2014/55), TS 1068	Gravimetrik	Kantitatif
119	Acılık Tayini	KAL/SOP85/PR15	Tahin ve Ürünleri	TS 2589	Organoleptik	Kalitatif
120	Alkolde Çözünmeyen Madde Tayini	KAL/SOP26/PR15	Konserve Gıdalar	TS ISO 23392, TS 382	Gravimetrik	Kantitatif
121	Asetil Metil Karbinol Tayini	KAL/SOP27/PR15	Sirke ve Ürünleri	TS 1880 EN 13188	Organoleptik	Kalitatif
122	Asit Sayısı ve Serbest Yağ Asitliği Tayini	KAL/SOP12/PR15	Hayvansal ve Bitkisel Katı ve Sıvı Yağlar (Ekstrakte Edilenler Dahil)	TS 894 TS EN ISO 660, TGK 2014/53	Titrimetrik	Kantitatif
123	Alkolde Çözünen Ekstrakt Tayini	KAL/SOP28/PR15	Baharatlar	TS 2135	Gravimetrik	Kantitatif

124	Karbondiyoksit Tayini	KAL/SOP29/PR15	Gazlı İçecekler, Alkollü İçkiler	TS 2259	Titrimetrik	Kantitatif
125	Çözünmeyen Safsızlıklar (Eterde Çözünmeyen Yabancı Madde) Tayini	KAL/SOP30/PR15	Hayvansal ve Bitkisel Yağlar	TS EN ISO 663	Gravimetrik	Kantitatif
126	Formol Sayısı Tayini	KAL/SOP31/PR15	Meyve Sebze suları	TS EN 1133	Titrimetrik	Kantitatif
127	Mineral Yağ Aranması	KAL/SOP33/PR15	Hayvansal ve Bitkisel Katı ve Sıvı Yağlar	TS 5039	Organoleptik	Kalitatif
128	Mineral Yağ Miktarı Tayini	KAL/SOP34/PR15	Hayvansal ve Bitkisel Katı ve Sıvı Yağlar	AOAC 966.17, TS 5039	Gravimetrik	Kantitatif
129	Nişasta Aranması	KAL/SOP20/PR15	Süt Tozları	EC 213-2001	Organoleptik	Kalitatif
130	Peroksidaz Testi	KAL/SOP35/PR15	Süt ve Süt Ürünleri, Dondurulmuş Bitkisel Ürünler	TS 1330, TS 1019, TS 7999	Organoleptik	Kalitatif
131	Peroksit Sayısı Tayini	KAL/SOP14/PR15	Hayvansal ve Bitkisel Katı ve Sıvı Yağlar (Ekstrakte Edilenler Dahil)	TS EN ISO 3960, TGK Zeytinyağı Ve Pirina Yağı Analiz Metotları Tebliği 2014/53, TS 2590, TS 894	Titrimetrik	Kantitatif
132	Toplam Şeker Tayini	KAL/SOP04/PR15	Gıdalar ve Yemler, Alkollü ve Alkolsüz İçecekler	TS 7780, TS 1466, EC 152/2009, TS 12232, TS 12001, TS 522, 21.01.2017 Tarih ve 29955 Sayılı Resmi Gazete 'de yayımlanan Yemlerin Resmî Kontrolü İçin Numune Alma ve Analiz Metotlarına Dair Yönetmelik	Titrimetrik	Kantitatif
133	Yüksek Alkol Uçucu Madde ve Metanol Miktarı Tayini	KAL/SOP40/PR15	Distile Alkollü ve Alkollü İçecekler	Türk Gıda Kodeksi Distile Alkollü İçkiler Analiz Metotları Tebliği (Tebliğ No: 2017/9)	Kromatografik	Kantitatif
134	Uçar Asit Tayini	KAL/SOP41/PR15	Şarap, Alkollü ve Alkolsüz İçecekler	TS 522, TS 6473	Titrimetrik	Kantitatif

135	Uçucu Olmayan Eter Ekstraktı Tayini	KAL/SOP42/PR15	Baharat ve Çeşni Veren Bitkiler, Renklendiriciler	TS 2137 ISO 1108	Gravimetrik	Kantitatif
136	Nişasta Tayini	KAL/SOP43/PR15	Gıdalar ve Yemler	TS 10520, RG 29955, EC 152/2009, TS ISO 6493	Polarimetrik	Kantitatif
137	Kolesterol Miktarı Tayini	KAL/SOP06/PR15	Gıdalar	J. Of AOAC Internatination, VOL 78, No:6,1995	Kromotografik	Kantitatif
138	Kokuşma Testi	KAL/SOP45/PR15	Et ve Et Ürünleri	TS 10580	Organoleptik	Kalitatif
139	Homojenizasyon Tayini	KAL/SOP46/PR15	Süt ve Süt Ürünleri	TS 1019	Gravimetrik	Kantitatif
140	Çözelti Rengi (ICUMSA) Tayini	KAL/SOP47/PR15	Şeker	TS 861	Spektrofotometrik	Kantitatif
141	Üre Tayini	KAL/SOP48/PR15	Yem ve Yem Katkı Maddeleri	RG 29955, TS 10338, AOAC 967.07	Spektrofotometrik	Kantitatif
142	Hazmolabilir Ham Protein Tayini	KAL/SOP49/PR15	Yem ve Yem Katkı Maddeleri	TS 12561, AOAC 971.09	Gravimetrik	Kantitatif
143	Üreaz Aktivitesi Tayini	KAL/SOP50/PR15	Soya ve Ürünleri	TS ISO 5506	Titrimetrik	Kantitatif
144	Amonyak Tayini	KAL/SOP53/PR15	Şarap	TS 522	Titrimetrik	Kantitatif
145	Asta Color Tayini	KAL/SOP55/PR15	Biber	AOAC 971.26	Spektrofotometrik	Kantitatif
146	Eterde Çözünmeyen Yabancı Maddeyle Kül Tayini	KAL/SOP58/PR15	Bitkisel Yağlar	TS 894	Gravimetrik	Kantitatif
147	Folin Ciocalteu İndeksi-Fenolik Maddeler Tayini	KAL/SOP60/PR15	Şarap	CONSLEG: 1990R2676-EN-Metot 41	Spektrofotometrik	Kantitatif
148	Furfural Testi	KAL/SOP61/PR15	Etil Alkol	CONSLEG-2000R1623-EN-Metot 11	Organoleptik	Kalitatif
149	Glikoz Şurubu Aranması	KAL/SOP62/PR15	Akide Şekeri	TS 7780	Organoleptik	Kalitatif

150	Kromatik Özellikler- Renk Yoğunluğu ve Renk Saflığı Tayini	KAL/SOP64/PR15	Şarap	CONSLEG:1990R 2676- Metot 40	Spektrofotometrik	Kantitatif
151	Pirina Yağı aranması	KAL/SOP66/PR15	Bitkisel ve Hayvansal Yağlar	TS 5042	Organoleptik	Kalitatif
152	Sedimentasyon Tayini	KAL/SOP67/PR15	Buğday, Buğday Unu	TS EN ISO 5529	Gravimetrik	Kantitatif
153	Serbest Amino Azot Tayini	KAL/SOP68/PR15	Bira	AOAC 977.08- AOAC 945.30	Spektrofotometrik	Kantitatif
154	Karbonhidrat Tayini	KAL/SOP70/PR15	Gıdalar	FAO Food and Nutrition Paper 77-Food Energy- Methods of Analysis and Conversion Factors	Gravimetrik	Kantitatif
155	Susam Yağı Aranması	KAL/SOP71/PR15	Hayvansal ve Bitkisel Yağlar	TS 5044	Organoleptik	Kalitatif
156	Hacmen Etil Alkol (Etanol) Analizi	KAL/SOP73/PR15	Distile Alkollü, Alkollü ve Alkolsüz İçecekler	CONSLEG 2000R2870, AOAC 984.14 AOAC 972.10 Türk Gıda Kodeksi Distile Alkollü İçkiler Analiz Metotları Tebliği (Tebliğ No: 2017/9)	Kromatografik	Kantitatif
157	Enerji Değeri ve Besin Öğesi Miktarlarının Hesaplanması Analizi	KAL/SOP74/PR15	Gıdalar	FAO Food and Nutrition Paper 77- Food Energy-Methods of Analysis and Conversion Factors, TGK Gıda Etiketleme ve Tüketicileri Bilgilendirme Yönetmeliği	Gravimetrik	Kantitatif
158	Metabolik Enerji Hesaplanması	KAL/SOP75/PR15	Yem ve Yem Katkı Maddeleri	TS 9610, EC 152/2009, RG 29955	Gravimetrik	Kantitatif
159	Siyanür Türevleri Tayini	KAL/SOP80/PR15	Şarap	CONSLEG: 1990R2676- EN-Metot 38	Titrimetrik	Kantitatif

160	p-Hidroksibenzoik Asit Esterleri Tayini (Metil p-hidroksi benzoat, Etil p-hidroksi benzoat, Propil p-hidroksi benzoat)	KKMT/SOP11/PR15	Tüm Gıdalar	NMKL 124.2	Kromatografik	Kantitatif
161	Organik Asit Tayini (Sitrik Asit, Askorbik Asit, Quinik Asit, Malik Asit, Laktik Asit, Propiyonik Asit, Asetik Asit, Suksinik Asit)	KKMT/SOP26/PR15	Alkolsüz İçecekler	İşletme İç Metot : "KKMT/SOP26/PR15" (AOAC 986.13)	Kromatografik	Kantitatif
162	Antioksidan Tayini (BHA, BHT, TBHQ, Propil Gallat)	KKMT/SOP28/PR15	Yağlar, Ekstrakte Edilmiş Yağlar	AOAC Official Method 983.15	Kromatografik	Kantitatif
163	Azodikarbonamit Analizi	KKMT/SOP29/PR15	Un ve Un Katkı Maddeleri	AACC Method 48-71A	Titrimetrik	Kantitatif
164	Potasyum hegzasiyanoferrat (II) (Potasyum ferrosiyanit) (Potasyum ferrosiyanür) Tayini	KKMT/SOP30/PR15	Tuz	EauSalt/AS 004-2008	Spekrofotometrik	Kantitatif
165	Formaldehit ve Hidrojen Peroksit Aranması	KKMT/SOP32/PR15	Çiğ Süt	TS 1018	Organoleptik	Kalitativ
166	Sülfat Tayini	KKMT/SOP33/PR15	Şarap	TS 522	Gravimetrik	Kantitatif
167	Karoten Tayini	KKMT/SOP34/PR15	Taze Meyve ve Sebzeler	İşletme İç Metot- "KKMT/SOP34/PR15" (AOAC Official Method 941.15)	Organoleptik	Kalitativ
168	Sukraloz Tayini	KKMT/SOP35/PR15	Gıda	TS EN 16155	Kromatografik	Kantitatif
169	Taurin Tayini	KKMT/SOP36/PR15	Enerji İçecekleri	Journal of Health Science, 54(6) 661-664 (2008)	Kromatografik	Kantitatif
170	Hidrojen Peroksit Tayini	KKMT/SOP37/PR15	Kuru İncir	İşletme İç Metot- "KKMT/SOP37/PR15" (TS 1018)	Organoleptik	Kalitativ

171	Karboksimetilselüloz Tayini	KKMT/SOP38/PR15	Süt ve Süt Ürünleri	İşletme İçi Metot- "KKMT/SOP38/PR15" (TS 6670)	Organoleptik	Kalitatif
172	Dithiocarbomate Grubu Pestisit Analizi	KKMT/SOP25/PR15	Sebze ve meyve (yüksek su içerikli ürünler, yüksek asit ve yüksek su içerikli ürünler), yüksek şeker ve düşük su içerikli ürünler, zor veya kendine özgü ürünler, yüksek yağ ve çok düşük su içerikli ürünler, yüksek yağ ve orta düzeyde su içerikli ürünler, yüksek nişasta ve/veya protein ve düşük su ve yağ içerikli ürünler	Acta Chim. Slov. 2006, 53, 100-104	Kromatografik	Kantitatif
173	DAS 40278-9 Mısır Tip Belirleme Analizi ' 2	MBAL/SOP58/PR15	Gıda ve Yem	EURL Metot (QT-EVE-ZM-004), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
174	Rutubet (Nem) Tayini – Kuru Madde Tayini	FAL/SOP41/ PR15	Mısır Cipsi ve Patates Cipsi, Kabak Çekirdeği	TS ISO 939	Volumetrik	Kantitatif
175	Rutubet (Nem) Tayini – Kuru Madde Tayini	FAL/SOP47/ PR15	Şeker ve Şekerli Ürünler, Meyve ve Sebze Ürünleri	TS 1129 ISO 1026, TGK İnsan Tüketimine Sunulan Şekerlerin Analiz Yöntemleri Tebliği	Gravimetrik	Kantitatif
176	Rutubet (Nem) Tayini – Kuru Madde Tayini	FAL/SOP03/ PR15	Baklagiller	TS ISO 24557	Gravimetrik	Kantitatif
177	Yağ Tayini (Gerber Metodu)	KAL/SOP08/PR15	Süt ve Süt Ürünleri	TS ISO 19662, TS 1018, TS 1330, TS ISO 3433, İşletme İçi Metot- KAL/SOP08/PR15 (Funke Gerber special method modifiye edilmiştir.)	Gravimetrik	Kantitatif
178	Yağ tayini	KAL/SOP87/PR15	Tereyağı	TS 1331	Gravimetrik	Kantitatif
179	Yağ tayini	KAL/SOP88/PR15	Bitkisel Margarin	TS 2812	Gravimetrik	Kantitatif

180	Diyet Lifi Tayini	KAL/SOP57/PR15	Gıdalar	Kit prosedürü, AOAC 991.43, AOAC 985.29 Sigma Kit Prosedürü	Gravimetrik	Kantitatif
181	Nem (Rutubet) Tayini – Kuru Madde Tayini	FAL/SOP01/PR15	Yem	RG 29955 Yemlerin Resmi Kontrolü İçin Numune Alma ve Analiz Metotlarına Dair Yönetmelik	Gravimetrik	Kantitatif
182	Organoleptik Muayene	FAL/SOP02/PR15	Tüm Gıdalar, Gıda Katkı Maddeleri, Yem ve Yem Katkı Maddeleri	TGK Tebliğ ve Yönetmelikleri, Ürün TS'leri, Food Chemical Codex, Ürün Spesifikasyonları, TS EN ISO 5492	Organoleptik	Kalitatif
183	Kutu Doldurma Oranı Tayini	FAL/SOP22/PR15	Taze Fasülye, Bezelye	TS 1468, TS 382	Gravimetrik	Kantitatif
184	Süzme Ağırlığı Tayini	FAL/SOP28/PR15	Taze Fasülye	TS 1468	Gravimetrik	Kantitatif
185	Yabancı Madde Tayini	FAL/SOP11/PR15	Tüm Gıdalar ve Yemler	TGK Tebliğleri, Ürün TS'leri İşletme İçi Metot “FAL/SOP11/PR15” (TGK Tebliğleri, Ürün TS'leri, Ürün Spektleri)	Gravimetrik	Kantitatif
186	Hacim Genmodsi Tayini	FAL/SOP44/PR15	Maraş Dondurması	TS 12649	Gravimetrik	Kantitatif
187	Suya Geçen Madde Miktarı Tayini	FAL/SOP25/PR15	Bulgur	TS 2284	Gravimetrik	Kantitatif
188	Sülfatlandırılmış Kül	FAL/SOP40/PR15	Nişasta ve Nişastalı Ürünler, Gıda ve Gıda Katkı Maddeleri, Şeker ve Şekerli Ürünler, Şekerlemeler, Şurup, Laktoz ve Dekstroz	TGK (Tebliğ no 2002/26), TS 6361 EN ISO 5809, Food Chemical Codex, Ürün Spesifikasyonları	Gravimetrik	Kantitatif
189	Yağ Asidi Metil Esterleri (Yağ Asitleri Kompozisyonu) Analizi	KAL/ SOP01/PR15	Hayvansal ve Bitkisel Katı ve Sıvı Yağlar (Ekstrakte Edilenler Dahil)	TS EN ISO 12966-1, TS EN ISO 12966-2, TS EN ISO 12966-4, TGK Zeytinyağı Ve Pirina Yağı	Kromatografik	Kantitatif

				Analiz Metotları Tebliği 2014/53		
C4:0 Methyl butyrate, C6:0 Methyl hexanoate, C8:0 Methyl octanoate, C10:0 Methyl decanoate, C11:0 Methyl undecanoate, C12:0 Methyl laurate, C13:0 Methyl tridecanoate, C14:0 Methyl myristate, C14:1 Methyl myristoleate (cis-9), C15:0 Methyl pentadecanoate, C15:1 Methyl pentadecanoate (cis-10), C16:0 Methyl palmitate, C16:1 Methyl palmitoleate (cis-9), C17:0 Methyl heptadecanoate, C17:1 Methyl heptadecenoate (cis-10), C18:0 Methyl stearate, C18:1 Methyl elaidate (trans-9), C18:1 Methyl oleate (cis-9), C18:2 Methyl linoleaidate (trans-9,12), C18:2 Methyl linoleate (cis-9,12), C20:0 Methyl arachidate, C18:3 Methyl γ -linolenate (cis-6, 9, 12), C20:1 Methyl eicosenoate (cis-11), C18:3 Methyl linolenate (cis-9,12,15), C21:0 Methyl heneicosanoate, C20:2 Methyl eicosadienoate (cis-11,14), C22:0 Methyl behenate, C20:3 Methyl eicosatrienoate (cis-8,11,14), C22:1 Methyl erucate (cis-13), C20:3 Methyl eicosatrienoate (cis-11,14,17), C20:4 Methyl arachidonate (cis-5,8,11,14), C23:0 Methyl tricosanoat, C22:2 Methyl docosadienoate (cis-13,16), C24:0 Methyl lignocerate, C20:5 Methyl eicosapentaenoate (cis-5,8,11,14,17), C24:1 Methyl nervonate (cis-15), C22:6 Methyl docosaheptaenoate (cis-4,7,10,13,16,19)						
190	Dekstroz Eşdeğeri Tayini	KAL/SOP72/PR15	Nişasta Hidroliz Ürünleri	TS EN ISO 5377, TGK 2002/26	Titrimetrik	Kantitatif
191	Doymuş Yağ Asitleri, Doymamış Yağ Asitleri, Trans Yağ Asitleri, Yağ Asidi Omega-3 (Alfa-Linolenik asit ALA, Eikosapentaenik asit EPA, Dokosaheptaenik asit DHA) Tayini	KAL/SOP86/PR15	Hayvansal ve Bitkisel Katı ve Sıvı Yağlar (Ekstrakte Edilenler Dahil)	TS EN ISO 12966-1, TS EN ISO 12966-2, TS EN ISO 12966-4 TGK Zeytinyağı ve Prina Yağı Analiz Metotları Tebliği (Tebliğ No: 2014/53)	Kromatografik	Kantitatif
192	Kolajen Analizi, Hidroksiprolin, Kolajen /Et Protein Oranı Tayini	KAL/SOP32/PR15	Et ve Et Ürünleri	TS 6236 ISO 3496, NMKL 127, TS 1748, ISO 937	Spektrofotometrik	Kantitatif
193	DAS-81419-2Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP66/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-014) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
194	DAS-44406-6 Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP67/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-015) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
195	DAS-68416-4 Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP71/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-013) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
196	Et Tür Tayini (Domuz DNA'sı Aranması) ¹⁻²	MBAL/SOP73/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Pork Assay)	Real Time PCR	Kalitatif
197	Et Tür Tayini (At DNA'sı Aranması) ¹⁻²	MBAL/SOP74/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Horse Assay)	Real Time PCR	Kalitatif

198	Et Tür Tayini (Eşek DNA'sı Aranması) ¹⁻²	MBAL/SOP75/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Donkey Assay)	Real Time PCR	Kalitatif
199	Et Tür Tayini (Sığır DNA'sı Aranması) ¹⁻²	MBAL/SOP76/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Beef Assay)	Real Time PCR	Kalitatif
200	Et Tür Tayini (Koyun DNA'sı Aranması) ¹⁻²	MBAL/SOP77/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Sheep Assay)	Real Time PCR	Kalitatif
201	Et Tür Tayini (Keçi DNA'sı Aranması) ¹⁻²	MBAL/SOP78/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Goat Assay)	Real Time PCR	Kalitatif
202	Et Tür Tayini (Tavuk DNA'sı Aranması) ¹⁻²	MBAL/SOP79/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Chicken Assay)	Real Time PCR	Kalitatif
203	Et Tür Tayini (Hindi DNA'sı Aranması) ¹⁻²	MBAL/SOP80/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Turkey Assay)	Real Time PCR	Kalitatif
204	Et Tür Tayini (Balık DNA'sı Aranması) ¹⁻²	MBAL/SOP81/PR15	Gıda, Yem, Jelatin ve Jelatin İçeren Ürünler	Kit Metodu (InnuDETECT Fish Assay)	Real Time PCR	Kalitatif
205	GDO Tarama Analizi (p35S, tNOS) ^{1 2}	MBAL/SOP82/PR15	Gıda ve Yem	EURL Metot (QL-ELE-00-013) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
206	GDO Tarama Analizi (pFMV) ^{1 2}	MBAL/SOP83/PR15	Gıda ve Yem	EURL Metot (QL-ELE-00-015) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
207	MON 87751 Soya Tip Belirleme Analizi ^{1 2}	MBAL/SOP86/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GM-016), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
208	Bitki Spesifik Kloroplast Geni Taraması ^{1 2}	MBAL/SOP87/PR15	Gıda ve Yem	İşletme İçi Metot (İnnuMix qPCR SyGreen Sensitive Modifiye Edilmiştir. MBAL/SOP87/PR15, ISO 21569, ISO 24276	Real Time PCR	Kalitatif
209	Bitki Spesifik Kloroplast Geni Taraması ^{1 2}	MBAL/SOP88/PR15	Gıda ve Yem	Kit Metodu (Bosphore Plant Detection Kit) ISO 21569, ISO 24276	Real Time PCR	Kalitatif
210	Rutubet (Nem) Tayini – Kuru Madde Tayini	FAL/SOP03/PR15	Hayvansal ve Bitkisel Katı ve Sıvı Yağlar, Meyve ve Sebze Suları, Çözünebilir Kahve	TS EN ISO 662, TS EN 12145, TS ISO 3726	Gravimetrik	Kantitatif

211	Rutubet (Nem) Tayini – Kuru Madde Tayini	FAL/SOP47/PR15	Tahin Helvası, Tahin, Lokum	TS 1201 EN ISO 1741, TS 2589, TS 2590, TS 8444/T1	Gravimetrik	Kantitatif
212	Bostwick Viskozitesi	FAL/SOP50/PR15	Salça, Ketçap, Sos ve Püreler, Meyve Şurubu vb. Ürünler	ASTM F-1080-93R02	Gravimetrik	Kantitatif
213	Gluten Analizi	MBAL/SOP85/PR15	Glutensiz ve Gluteni Azaltılmış Tüm Gıdalar	Kit Metodu (AgraQuant Gluten G12), AOAC 2014.03	ELISA	Kantitatif
214	Polar Madde Tayini	KAL/SOP92/PR15	Hayvansal ve Bitkisel Yağlar, Ekstrakte Edilen Yağlar	TS EN ISO 8420	Gravimetrik	Kantitatif
215	Prolin Tayini	KAL/SOP91/PR15	Meyve Sebze Suları	TS EN 1141	Spektrofotometrik	Kantitatif
216	Sterol, Eritrodiol, Uvaol Bileşimi Ve Miktarı Tayini	KAL/SOP105/PR15	Bitkisel Ve Hayvansal Yağlar, Ekstrakte Edilen Yağlar	TGK Zeytinyağı Ve Pirina Yağı Analiz Metotları Tebliği 2014/53, ISO 12228-1, ISO 12228-2	Kromatografik	Kantitatif
Kolesterol, Brassikasterol, Ergosterol, 24-metilen-kolesterol, Kampesterol, Kampestanol, Stigmasterol, Delta-7-kampesterol, Delta-5,23-stigmastadienol, Klerosterol, Beta-sitosterol, Sitostanol, Delta-5-avenasterol, Delta-5.24-stigmastadienol, Delta-7-stigmastenol, Delta-7-avenasterol, Eritrodiol, Uvaol, Toplam Sterol, Toplam Eritrodiol + Uvaol, Toplam Beta-Sitosteoller						
217	UV Özgül Absorbans Değeri	KAL/SOP94/PR15	Zeytinyağı-Pirinayağı	TGK Zeytinyağı Ve Pirina Yağı Analiz Metotları Tebliği 2014/53- COI T.20 doc 19- TS EN ISO 3656 A1- TS 342	Spektrofotometre	Kantitatif
218	Susam Yağı Miktarı	KAL/SOP93/PR15	Tahin Helvası, Tahin	TS 2590, TS 2589	Gravimetrik	Kantitatif
219	Süt Yağı Harici Yağ Aranması	KAL/SOP90/PR15	Süt Ve Süt Ürünleri (olgunlaştırılmış peynir ve taze peynir hariç)	TS EN ISO 17678	Kromatografik	Kantitatif
220	Bitkisel Yağlarda Toplam Migrasyon Analizi Metot 1- Tam Daldırma Metodu	KAL/SOP113/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-2 Metot 1	Kromatografik	Kantitatif
221	Bitkisel Yağlarda Toplam Migrasyon	KAL/SOP113/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-2	Kromatografik	Kantitatif

	Analizi Metot 2- Hücre Metodu			Metot 2		
222	Bitkisel Yağlarda Toplam Migrasyon Analizi Metot 3- Torba Metodu	KAL/SOP113/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-2 Metot 3	Kromotografik	Kantitatif
223	Bitkisel Yağlarda Toplam Migrasyon Analizi Metot 4- Ters Torba Metodu	KAL/SOP113/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-2 Metot 4	Kromotografik	Kantitatif
224	Bitkisel Yağlarda Toplam Migrasyon Analizi Metot 5- Doldurma Metodu	KAL/SOP113/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-2 Metot 5	Kromotografik	Kantitatif
225	Buharlaşılabilen Gıda Benzerlerinde Toplam Migrasyon Analizi Metot 1- Tam Daldırma Metodu	KAL/SOP114/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-3 Metot 1	Gravimetrik	Kantitatif
226	Buharlaşılabilen Gıda Benzerlerinde Toplam Migrasyon Analizi Metot 2- Hücre Metodu	KAL/SOP114/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-3 Metot 2	Gravimetrik	Kantitatif
227	Buharlaşılabilen Gıda Benzerlerinde Toplam Migrasyon Analizi Metot 3- Torba Metodu	KAL/SOP114/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-3 Metot 3	Gravimetrik	Kantitatif
228	Buharlaşılabilen Gıda Benzerlerinde Toplam Migrasyon Analizi Metot 4- Ters Torba Metodu	KAL/SOP114/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-3 Metot 4	Gravimetrik	Kantitatif
229	Buharlaşılabilen Gıda Benzerlerinde Toplam Migrasyon Analizi Metot 5- Doldurma Metodu	KAL/SOP114/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-3 Metot 5	Gravimetrik	Kantitatif
230	Toplam Migrasyon Analizi Yüksek Sıcaklık Yerine Geçme Metodu (Yağlı Gıdalar) Modifiye	KAL/SOP100/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, TS EN 1186-13	Gravimetrik	Kantitatif

	Polifenilen Oksit (MPPO)					
231	Toplam Migrasyon Analizi Polifenilen Oksit (MPPO)	KAL/SOP100/PR15	Gıda ile Temas Eden Kâğıt ve Karton Malzemeler	TS EN 14338	Gravimetrik	Kantitatif
232	Yapı Tayini FT-IR Metodu	KAL/SOP107/PR15	Gıda ile Temas Eden Madde ve Malzemeler	İşletme İçi Metot "KAL/SOP107/PR15" (Shimadzu Library Book)	FT-IR/ATR	Kalitativ
233	Bacillus cereus Sayımı	MAL/SOP58/PR15	Tüm Gıda ve Yem	COMPASS Bacillus Reference BKR 23/06-02/10 Detection and enumeration of Bacillus cereus group in human and animal food product.	Yayma / Dökme Plak	Kantitatif
234	Koagülaz pozitif Staphylococci (Staphylococcus aureus ve diğer türleri) Sayımı	MAL/SOP59/PR15	Tüm Gıda	Easy Staph Reference BKR 23/10-12/15 Enumeration Of Coagulase Positive Staphylococci group in human food product.	Yayma / Dökme Plak	Kantitatif
235	Salmonella spp. Aranması	MAL/SOP60/PR15	Tüm Gıda ve Yem	IRIS Salmonella Reference BKR 23/07-10/11 Detection of Salmonella spp. group in human and animal food product.	Yayma / Dökme Plak	Kalitativ
236	Listeria mono Aranması	MAL/SOP63/PR15	Tüm Gıda	Compass Listeria Reference BKR 23/02-11/02 Detection of Listeria monocytogenes group in human food product.	Yayma / Dökme Plak	Kalitativ

237	Maya Küf Sayımı	MAL/SOP61/PR15	Tüm Gıda ve Yem	SYMPHONY Agar reference BKR 23/11-12/18(Enumeration of yeast and molds in all human and animal food products)	Yayma / Dökme Plak	Kantitatif
238	Stafilokokkal Enterotoksin Aranması	MAL/SOP62/PR15	Tüm Gıda	EU-RL European Screening Method Ver.5 Ridascreeen Set Total Kit R 4105	Kit Metodu	Kalitatif
239	Baryum (Ba), Kobalt (Co), Bakır (Cu), Demir (Fe), Lityum (Li), Mangan (Mn), Çinko (Zn), Alüminyum (Al), Nikel (Ni) Analizi	KKMT/SOP47/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler	TS EN 1186-1, ISO 17294-1, ISO 17294-2, TS EN 13130-1	ICP-MS	Kantitatif
240	Arsenik (As), Cıva (Hg), Kadmiyum (Cd), Kurşun (Pb) Kalay (Sn) Sodyum (Na), Çinko (Zn), Kalsiyum (Ca), Magnezyum (Mg), Potasyum (K) Bakır (Cu), Krom (Cr) Demir (Fe), Nikel (Ni), Manganez (Mn), Fosfor (P) ve Selenyum Analizi	KKMT/SOP13-2/PR15	Tüm Gıdalar (Şarap Hariç)	NMKL 161, NMKL 186, NKML 191, TS EN 15621, TS EN 15510	ICP-MS	Kantitatif

241	Arsenik (As), Cıva (Hg), Kadmiyum (Cd), Kurşun (Pb) Kalay (Sn) Sodyum (Na), Çinko (Zn), Kalsiyum (Ca), Magnezyum (Mg), Potasyum (K) Bakır (Cu), Krom (Cr) Demir (Fe), Nikel (Ni), Manganez (Mn), Fosfor (P) ve Selenyum Analizi	KKMT/SOP13-2/PR15	Tüm Gıda Katkıları	NMKL 161, NMKL 186, NKML 191, TS EN 15621, TS EN 15510	ICP-MS	Kantitatif
242	Arsenik (As), Cıva (Hg), Kadmiyum (Cd), Kurşun (Pb) Kalay (Sn) Sodyum (Na), Çinko (Zn), Kalsiyum (Ca), Magnezyum (Mg), Potasyum (K) Bakır (Cu), Krom (Cr) Demir (Fe), Nikel (Ni), Manganez (Mn), Fosfor (P) ve Selenyum Analizi	KKMT/SOP13-2/PR15	Tüm Yemler	NMKL 161, NMKL 186, NKML 191, TS EN 15621, TS EN 15510	ICP-MS	Kantitatif
243	Arsenik (As), Cıva (Hg), Kadmiyum (Cd), Kurşun (Pb) Kalay (Sn) Sodyum (Na), Çinko (Zn), Kalsiyum (Ca), Magnezyum (Mg), Potasyum (K) Bakır (Cu), Krom (Cr) Demir (Fe), Nikel (Ni), Manganez	KKMT/SOP13-2/PR15	Tüm Yem Katkıları	NMKL 161, NMKL 186, NKML 191, TS EN 15621, TS EN 15510	ICP-MS	Kantitatif

	(Mn), Fosfor (P) ve Selenyum Analizi					
244	Spesifik Migrasyon Primer Aromatik Aminlerin (PAA) Analizi	KKMT/SOP45/PR15	Gıda ile Temas Eden Madde ve Malzemeler	EUR 24815, TS EN 13130-1	Kromatografik	Kantitatif
2,4,5-Trimethylaniline (2,4,5-TMA), 2,4-Dimethylaniline(2,4-DMA), 2,4-Toluenediamine (2,4-TDA), 2,6-Dimethylaniline (2,6-DMA), 2,6-Toluenediamine (2,6-TDA), 2-Methoxy-5-Methylaniline (2-M-5-MA), 3,3-Dimethylbenzidine (3,3-DMB), 4,4'-Methylenedianiline (4,4'-MDA), 4,4-Diaminodiphenylether (4,4'-DPE), 4-Chloro-Aniline (4-CA),4-Chloro-o-Toluidine (4-CoT),4-Methoxy-m-phenylenediamine (4-M-mPDA), Aniline(ANL), Benzidine (BNZ), o-Anisidine (o-ASD), o-Toluidine (o-T)						
245	Spesifik Migrasyon Fitalat Migrasyon Tayini	KKMT/SOP43/PR15	Gıda ile Temas Eden Madde ve Malzemeler	İşletme içi Metot (TS EN 14372'den modifiye), TS EN 13130-1, Food Additives and Contaminants, Report of an interlaboratory comparison organised by the EU Reference Laboratory for Food Contact Materials. (EUR 24387 EN 2010),• Türk Gıda Kodeksi Gıda İle Temas Eden Plastik Madde Ve Malzemeler Tebliği. (Tebliğ No: 2019/44)	Kromatografik	Kantitatif
Phthalic acid, bis-iso-decyl ester(DIDP),Phthalic acid, benzylbutyl ester(BBP), phthalic acid, bis-iso-butyl ester (DIBP), phthalic acid, bis-butyl ester (DBP), Bis (2-ethylhexyl) phthalate(DEHP), Diisononyl phthalate(DINP)						

246	Spesifik Migrasyon (2,2-bis(4-hydroxyphenyl) propane (Bisphenol A (BPA)) Analizi	KKMT/SOP44/PR15	Gıda ile Temas Eden Plastik Madde ve Malzemeler (Çocuklar için Kullanım ve Bakım Eşyaları, Plastik Çatal, Bıçak ve Mama Eşyaları (biberon vb.) ve Polikarbonat Türü Plastikler Kapsamına Giren Malzemeler)	EN 14372, 5991-5037 EN	Kromatografik	Kantitatif
247	Akrilamid Analizi	KKMT/SOP41/PR15	Gıda, Bebek Maması, Bebek ve Küçük Çocuk Ek Gıdası, Devam Formülleri	İşletme İçi Metot Agilent (Quantification of Acrylamide in a Variety of Food Matrices by LC/MS/MS Triple Quadrupole), Akrilamid Yönetmeliği	Kromatografik	Kantitatif
248	Kinin Analizi	KKMT/SOP42/PR15	Alkolsüz İçecekler, Enerji ve Sporcu İçecekleri	TS 12625	Kromatografik	Kantitatif
249	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Bifenthrin, Biphenyl, Bixafen, Bromocyclen, Bromophos Ethyl, Bromopropylate, Captan, Carbophenothion, Chinomethionate, Chlorbenseide, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenprop Methyl, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Cypermethrin-Alpha, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dichlorprop-ethyl; Dichlorprop-methyl; Dicloran, Dicofof, Dieldrin, Dimethipin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etridiazole, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoprene, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Captan met. THPI; Trifluralin, Vinclozolin						

250	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Asit Ve Yüksek Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Bifenthrin, Biphenyl, Bixafen, Bromocyclen, Bromophos Ethyl, Bromopropylate, Captan, Carbophenothion, Chinomethionate, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenprop Methyl, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Cypermethrin-Alpha, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dichlorprop-ethyl; Dichlorprop-methyl; Dicloran, Dicofo, Dieldrin, Dimethipin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etridiazole, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoprene, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Captan met. THPI; Trifluralin, Vinclozolin						
251	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Şeker Ve Düşük Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
252	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Yağ Ve Çok Düşük Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						

Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrpyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
253	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Yağ Ve Orta Düzeyde Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrpyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
254	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Nişasta Ve/Veya Protein Ve Düşük Su Ve Yağ İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrpyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
255	Pestisit Analizi	KKMT/SOP01/PR15	Zor Veya Özgün Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin						

aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrpyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
256	Pestisit Analizi	KKMT/SOP01/PR15	Et (Kas) Ve Su Ürünleri	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrpyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
257	Pestisit Analizi	KKMT/SOP01/PR15	Hayvansal Kökenli Gıdalardan Elde Edilen Yağlar	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrpyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
258	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Su İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin,						

Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofol, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
259	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Asit Ve Yüksek Su İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofol, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
260	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Yağ/ Katı Yağ Ve Düşük Su İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofol, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
261	Pestisit Analizi	KKMT/SOP01/PR15	Orta Yağ Ve Düşük Su İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone,						

Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrifos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
262	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Nişasta Ve/Veya Protein Ve Düşük Su Ve Yağ İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrifos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
263	Pestisit Analizi	KKMT/SOP01/PR15	Zor Veya Özgün Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS : 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrifos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
264	Pestisit Analizi	KKMT/SOP01/PR15	Et Ve Su Ürünleri (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif

GC-MS/MS: 1-Naphtylacetamide, 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
265	Pestisit Analizi	KKMT/SOP01/PR15	Süt Ve Süt Ürünleri (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrapyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						
266	Pestisit Analizi	KKMT/SOP21/PR15	Bebek Formülleri ve Devam Formülleri, Bebek ve Küçük Çocuk Ek Gıdaları	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 2,4-DDD, 2,4-DDE, 2,4-DDT, 2-Chloroaniline, 2-Phenylphenol, 3-Chloroaniline, 4,4-DDD, 4,4-DDE, 4,4-DDT, 4-Chloroaniline, Aclonifen, Acrinathrin, Aldrin, Anilazine, Benfluralin, Biphenyl, Bromocyclen, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlordane-cis, Chlordane-trans, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofo, Dieldrin, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos,						

Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrpyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin

267	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
-----	------------------	-----------------	----------------------------	--	---------------	------------

LC-MS/MS: 2,4 DB; 2,4 DB methylester; 2,4-5-T; 2,4-D; 2-hydroxy propycarbazone; 3,4,5-Trimethacarb; 4 nitrophenol sodium; 4-(3-indolyl)-butric acid (Indolylbutyric acid); 4-Aminopyridine; 5 nitroguaiacol sodium salt; 6-benzylaminopurine (6-benzyladenine); 6-chloronicotinic acid; Avermectin B1a; Avermectin B1b; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Allidochlor; Ametoctradin; Ametryn; Amidosulfuron; Aminocarb; Aminopyralid; Amisulbrom; Amitraz; Ancymidol; Anilofos; Aramite; Asulam; Atraton; Atrazine; Azadirachtin; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoconazole; Azoxystrobin; Barban; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bensulide; Bentazone; Benthiavalicarb-Isopropyl; Benzobicyclon; Benzoximate; Bifenazate; Bifenox; Bispyribac; Bispyribac Sodium; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butachlor; Butafenacil; Butoxycarboxim; Butralin; Buturon; Butylate; Cadusafos; Carbaryl; Carbazim+Benomyl; Carbetamide; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Carpromaid; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorflazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chlorthiophos; Chlortoluron; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Clomazone; Cloquintocet-1-methylhexylester; Cloransulam-methyl; Clothianidin; Coumaphos; Cruformate; Cyanazine; Cyantraniliprole; Cyazofamid; Cyclanilide; Cycloate; Cyclosulfamuron; Cycloxydim; Cycluron; Cyflufenamid; Cyflumetofen; Cymoxanil; Cyprazine; Cyproconazole; Cyprodinil; Cyprosulfamide; Cyromazine; Daminozide; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Desmetryn+Simetryn; Diafenthiuron; Dialifos; Diallyate; Diamuron; Diazinon; Dicamba; Dichlofluanid; Dichlormid; Dichlorprop-P+ Dichlorprop; Dichlorvos; Diclobutrazol; Diclofop-Methyl; Diclosulam; Dicrotophos; Dicyclanil; Diethofencarb; Difenconazole; Difenoxuron; Diflubenzuron; Diflufenican; Dimefox; Dimefuron; Dimethacarb; Dimethachlor; Dimethametryn; Dimethenamid; Dimethilan; Dimethoate; Dimethomorph; Dimethomorph E; Dimethomorph Z; Dimoxystrobin; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Disulfoton sulphone; Disulfoton sulfoxide; Ditalimfos; Dithianon; Dithiopyr; Diuron; Amitraz met. DMF; DNOC; Dodemorph; Dodine; Edifenphos; Ehtiofencarb sulphone; EmamectinB1A; EPN; Epoxiconazole; EPTC; Esprocarb; Etaconazole; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etobenzanid; Etofenprox; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenfuram; Fenhexamid; Fenobucarb; Fenoxanil; Fenoxaprop+Fenoxaprop-P; Fenoxaprop-ethyl+Fenoxaprop-P-ethyl; Fenoxycarb; Fenpiclonil; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyrazimine; Fenpyroximate E; Fensulfothion; Fensulfothion oxon; Fensulfothion sulfone; Fenthion; Fenthion oxon; Fenthion sulfone; Fenthion sulfoxide; Fentin+Fentin Chloride; Fenuron; Ferimzone; Fipronil; Fipronil desulfenil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Flonicamid; Florasulam; Florryauxifen- benzyl; Fluacrypyrim; Fluazifop+Fluazifop-P; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubendiamide; Flubenzimine; Flucarbazon; Flucetosulfuron; Flucycloxuron; Flucythrinate; Fludioxonil; Flufenacet; Flufenacet alcohol; Flufenacet thioglycolate sulfoxide; Flufenacet oxalate; Flufenoxuron; Flufenzine; Flumetsulam; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Flupyradifurone; Fluquinconazole; Fluridone; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Fluthiacet methyl; Flutolanil; Flutriafol; Fluxapyroxad; Fomesafen; Fonofos; Foramsufuron; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furalaxyl; Furametpyr; Furathiocarb; Furmecyclox; Halauxifen-methyl; Halosulfuron methyl; Haloxifop; Haloxifop-2-Ethoxy-Ethyl(etotyl); Haloxifop-methyl; Haloxifop-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexazinone; Hexythiazox; Hydroxyquinoline; Imazalil; Imazamethabenz-methyl; Imazamox; Imazapic; Imazapyr; Imazaquin; Imazethapyr; Imazosulfuron; Imibenconazole; Imidacloprid; Inabenfide; Indanofan; Indaziflam; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprobenphos; Iprodione; Iprovalicarb; Isazofos; Isocarbamide; Isocarbophos; Isofenphos; Isofenphos-methyl; Isonoruron; Isoprocab; Isopropalin; Isoprothiolane; Isoproturon; Isopyrazam; Isoxaben; Isoxaflutole; Isoxathion; Karbutilate; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mefenacet; Mefentrifluconazole; Mepanipyrim; Mephosfolan; Mepronil; Mesosulfuron-methyl; Mesotrione; Metaflumizone; Metalaxyl+Metalaxyl-M; Metamitron; Metazachlor; Metconazole; Metconazole Cis; Methabenzthiazuron; Methacifos; Methamidophos; Methamitron; Methfuroxam; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methopyrotryne; Methosfolan; Methoxyfenozone; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Mexacarbate; Milbectin A3; Milbectin A4; Molinate; Monocrotophos; Monolinuron;

Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Napropamide; Naptalam; Neburon; Nicosulfuron; Nitenpyram; Nitratin; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orbencarb; Orthosulfamuron; Oryzalin; Oxadiargyl; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Oxydemeton-methyl; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pencycuron-PB-amine; Pendimethalin; Penflufen; Penoxsulam; Penthioopyrad; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picloram; Picolinafen; Picoxystrobin; Pinoxaden; Piperonyl butoxide; Piperophos; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Pretilachlor; Primisulfuron methyl; Prochloraz; Profenofos; Profoxydim; Prohexadione; Promecarb; Prometon; Prometryn; Propachlor; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propoxur; Propoxycarbazon; Propyzamide; Proquinazid; Prosulfocarb; Prosulfuron; Prothioconazole; Prothioconazole desthio; Prothiofos; Pymetrozine; Pyracarbolid; Pyraclofos; Pyraclostrobin; Pyraflufen ethyl; Pyraflufen monohydrate; Pyrasulfotole; Pyrazophos; Pyrethrins (Pyrethrin I+II); Pyributicarb; Pyridaben; Pyridalyl; Pyridaphenthion; Pyridate; Pyrifeno; Pyrimethanil; Pyrimidifen; Pyriofenone; Pyriproxyfen; Pyriothibac; Pyroquilon; Pyroxasulfone; Pyroxsulam; Quinalphos; Quinclorac; Quinmerac; Quinoclamine; Quinoxifen; Quizalofop; Quizalofop-methyl; Quizalofop-P-ethyl; Quizalofop-P-tefuryl; Rimsulfuron; Rotenone; Sebutylazir; Sedaxane; Sethoxydim Solution; Siduron; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spiroclufen; Spirotetramat; Spirotetramat-enol; Spiroxamine; Sulcotrione; Sulfosulfuron; Sulfotep; Sulfoxaflor; Sulprofos; Tau-fluvalinate; Tebuconazole; Tebufenozide; Tebufenpyrad; Tebufenozide; Tebupirimfos; Tebutam; Tebutiuron; Teflubenzuron; Tembotrione; Temephos; Tepraloxym; Tepraloxym 3-hydroxy-glutaric acid; Terbufos sulfone; Terbufos sulfoxide; Terbumeton; Terbumeton desethyl; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiacloprid; Thiamethoxam; Thiazopyr; Thidiazuron; Thiencazone methyl; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thionazin; Thiophanate; Thiophanate-Methyl; Tiadinil; Tolclofos-Methyl; Tolfenpyrad; Tolyfluand; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazamate; Triazophos; Tribenuron-methyl; Trichlorfon; Triclopyr; Tricyclazole; Tridemorph; Trietazine; Trifloxystrobin; Trifloxysulfuron; Triflumizole; Triflumizole metabolite FM-6-1; Triflururon; Triflursulfuron; Triflursulfuron-methyl; Triforine; Trinexapac methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Valifenalate; Vamidothion; Vernolat; Zoxamide;

268	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Asit ve Yüksek Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
-----	------------------	-----------------	---	---	---------------	------------

LC-MS/MS: 2,4 DB; 2,4 DB methylester; 2,4-5-T; 2,4-D; 2-hydroxy propycarbazone; 3,4,5-Trimethacarb; 4 nitrophenol sodium; 4-(3-indolyl)-butric acid (Indolylbutyric acid); 4-Aminopyridine; 5 nitroguaiacol sodium salt; 6-benzylaminopurine (6-benzyladenine); 6-chloronicotinic acid; Avermectin B1a; Avermectin B1b; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Allidochlor; Ametoctadin; Ametryn; Amidosulfuron; Aminocarb; Aminopyralid; Amisulbrom; Amitraz; Ancymidol; Anilofos; Aramite; Asulam; Atraton; Atrazine; Azadirachtin; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoconazole; Azoxystrobin; Barban; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bensulide; Bentazone; Benthiavalicarb-Isopropyl; Benzobicyclon; Benzoximate; Bifenazate; Bifenox; Bifenthrin; Bispyribac; Bispyribac Sodium; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butachlor; Butafenacil; Butoxycarboxim; Butralin; Buturon; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbetamide; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Carpomaid; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chlorthiophos; Chlortoluron; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Clomazone; Cloquintocet-1-methylhexylester; Cloransulam-methyl; Clothianidin; Coumaphos; Cruformate; Cyanazine; Cyantraniliprole; Cyazofamid; Cyclanilide; Cycloate; Cyclosulfamuron; Cycloxydim; Cycluron; Cyflufenamid; Cyflumetofen; Cymoxanil; Cyprazine; Cyproconazole; Cyprodinil; Cyprosulfamide; Cyromazine; Daminozide; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Desmetryn+Simetryn; Diafenthiuron; Dialifos; Diallylate; Diamuron; Diazinon; Dicamba; Dichlofluand; Dichlormid; Dichlorprop-P+ Dichlorprop; Dichlorvos; Diclobutrazol; Diclofop-Methyl; Diclosulam; Dicrotophos; Dicyclanil; Diethofencarb; Difenconazole; Difenoxuron; Diflubenzuron; Diflufenican; Dimefox; Dimefuron; Dimethacarb; Dimethachlor; Dimethametryn; Dimethenamid; Dimethilan; Dimethoate; Dimethomorph; Dimethomorph E; Dimethomorph Z; Dimoxystrobin; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Disulfoton sulphone; Disulfoton sulfoxide; Ditalimfos; Dithianon; Dithiopyr; Diuron; Amitraz met. DMF; DNOC; Dodemorph; Dodine; Edifenphos; Ehtiofencarb sulphone; EmamectinB1A; EPN; Epoxiconazole; EPTC; Esprocarb; Etaconazole; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etobenzanid; Etofenprox; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenfuram; Fenhexamid; Fenobucarb; Fenoxanil; Fenoxaprop+Fenoxaprop-P; Fenoxaprop-ethyl+Fenoxaprop-P-ethyl; Fenoxycarb; Fenpiclonil; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproxymate;

Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorflazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; -S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallylate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenoconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafole; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetenate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxyfop; Haloxyfop-2-Ethoxy-Ethyl(etotyl); Haloxyfop-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarboxiphos; Isopropyl; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metconazole; Methacryfos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxyfen; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepaloxymid; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiadiazuron; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflurosulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide

270	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Yağ Ve Çok Düşük Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
-----	------------------	-----------------	---	---	---------------	------------

LC-MS/MS: 2,4-D; 2,4-5-T; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthion; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobil; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetenate ;Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M;Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraon-Ethyl; Paraon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiacloprid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflurosulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide;

271	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Yağ Ve Orta Düzeyde Su İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
-----	------------------	-----------------	--	---	---------------	------------

LC-MS/MS; 2,4-D; 2,4-5-T; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorflazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthion; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobil; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetenate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isafenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraon-Ethyl; Paraon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxymid; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiacloprid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflusulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide;

272	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Nişasta Ve/Veya Protein Ve Düşük Su Ve Yağ İçerikli Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-D; 2,4-5-T; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Bflubutamid; Benalaxyl; Bendiocarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoxazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetenate; Fosthiazate; Fuberidazole; Halosulfuron methyl; Haloxyfop; Haloxyfop-2-Ethoxy-Ethyl(etotyl); Haloxyfop-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam;—Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuairimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Promeytryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxyfen; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiacloprid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-						

Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflururon; Triflurosulfuron-methyl; Trinexapac-ethyl; Triticonazole; Vamidothion; Zoxamide;						
273	Pestisit Analizi	KKMT/SOP01/PR15	Zor Veya Özgün Ürünler	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-D; 2,4-5-T; 3,4,5-Trimethacarb; Acephate; Acetamid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Beflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthialacarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethamsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobil; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isoprazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metconazole; Methacifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuairimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Prophanil; Propiconazole; Propoxur; Propylamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quiazaalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufenozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiacloprid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-						

sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflusulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide;						
274	Pestisit Analizi	KKMT/SOP01/PR15	Et (Kas) Ve Su Ürünleri	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-D; 2,4-5-T; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Beflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthion; Dialifos; Diallylate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobil; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; -Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quiazaflufop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxyfen; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos;						

Tetraconazole; Tetramethrin; Thiabendazole; Thiacloprid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflusulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide;						
275	Pestisit Analizi	KKMT/SOP01/PR15	Süt ve Süt Ürünleri	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Beflubutamid; Benalaxyl; Bendiocarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproxymate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Fluquinconazole; Flurochloridone; Fluroxyppy; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Halosulfuron methyl; Haloxyfop; Haloxyfop-2-Ethoxy-Ethyl(etotyl); Haloxyfop-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaaxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxyfen; Quizalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiacloprid; Thiamethoxam;						

Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluand; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflusulfuron-methyl; Trinexapac-ethyl; Triticonazole; Vamidothion; Zoxamide;						
276	Pestisit Analizi	KKMT/SOP01/PR15	Hayvansal Kökenli Gıdalar Elde Edilen Yağlar	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Beflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthialdicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorflazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isofenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanol; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Prophanil; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quazalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Quizalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufenozide; Teflubenzuron;						

Tepaloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiacloprid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluamid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflusulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide;						
277	Pestisit Analizi	KKMT/SOP21/PR15	Bebek Formülleri ve Devam Formülleri, Bebek ve Küçük Çocuk Ek Gıdaları	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4,5-T; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-methyl; Alachlor; Aldicarb; Aldicarb-sulfone; Aldicarb-sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-ethyl; Azinphos-methyl; Aziprotryne; Azoxystrobin; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Benomyl; Bensulfuron-methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim-sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-methyl sulfone; Demeton-S-methyl sulfoxide (Oxydemeton-methyl); Desmedipham; Diafenthiuron; Dialifos; Diallylate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Disulfoton sulfone; Disulfoton sulfoxide; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoxazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos-sulfone; Fenamiphos-sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin-oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenpyroximate-(E); Fensulfothion; Fensulfothion-oxon; Fensulfothion-sulfone; Fenthion; Fenthion-sulfone; Fenthion-sulfoxide; Fentin oxid; Fentin+Fentin Chloride; Fipronil; Fipronil-desulfinyl; Fipronil-sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetenate(hydrochloride); Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron-methyl; Haloxyfop+Haloxypyr-R; Haloxypyr+Haloxypyr-R; Haloxypyr-methyl+Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-methyl; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarboxiphos; Isopropyl; Isoproturon; Isopropyl; Isoxaflutole; Kresoxim-methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metconazole; Methacryfos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb-sulfone; Methiocarb-sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; N-(2,4-Dimethylphenyl)Formamide; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace;						

Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate-sulfone; Phorate-sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propoxur; Propylenethiourea (PTU); Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionat (Chinomethionat); Quinoxifen; Rimsulfuron; Sethoxydim; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-Fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbufos-Sulfone; Terbufos-Sulfoxide; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiachlorid; Thiamethoxam; Thidiazuron; Thifensulfuron-Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-methyl; Tolclofos-methyl; Tolfenpyrad; Tolyfluanid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflururon; Triflurosulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide;

278	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Su İçerikli Ürünler(Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
-----	------------------	-----------------	---------------------------------	---	---------------	------------

LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Beflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorflazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthion; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethamsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoxazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Fibronil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobil; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isoprazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron;

Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pymethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Quizalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiachlopid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflusaluron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide						
279	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Asit Ve Yüksek Su İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Befulbutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluaazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthion; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Fibronil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxyfop; Haloxyfop-2-Ethoxy-Ethyl(etotyl); Haloxyfop-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron;						

Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; ; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxsulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Quizalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufenozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiachlorid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflururon; Triflursulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide						
280	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Yağ/ Katı Yağ Ve Düşük Su İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Beflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Fibronil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastobin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isoprazam; Isoxaflutole; Kresoxim-Methyl; Lenacil;						

Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; ; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Quizalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiachlorid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluamid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflururon; Triflursulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide						
281	Pestisit Analizi	KKMT/SOP01/PR15	Orta Yağ Ve Düşük Su İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Fibronil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Fluquinconazole; Flurochloridone; Fluroxyppyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxyfop; Haloxyfop-2-Ethoxy-Ethyl(etotyl); Haloxyfop-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid;						

Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isafenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaaxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; ; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuairimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Prophan; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quazalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxyfen; Quizalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufenozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiachlorid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluand; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflururon; Triflursulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide						
282	Pestisit Analizi	KKMT/SOP01/PR15	Yüksek Nişasta Ve/Veya Protein Ve Düşük Su Ve Yağ İçerikli Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Pestisit Analizi
LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamidrid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorprophan; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluand; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Fibronil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole;						

Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxyfop; Haloxyfop-2-Ethoxy-Ethyl(etotyl); Haloxyfop-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isopyrazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; -Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Propham; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quaizalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Quizalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiachloprid; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflumuron; Triflurosulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide						
283	Pestisit Analizi	KKMT/SOP01/PR15	Zor Veya Özgün Ürünler (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Pestisit Analizi
LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Fibronil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil;						

Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobilin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Iproconazole; Iprodione; Iprovalicarb; Isocarbophos; Isafenphos; Isoproturon; Isoprazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Meposfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacryfos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuaimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Prophan; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quiazaflor-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Quiazaflor-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufenozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiadiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluamid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflururon; Triflurosulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide

284	Pestisit Analizi	KKMT/SOP01/PR15	Et Ve Su Ürünleri (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
-----	------------------	-----------------	-------------------------	---	---------------	------------

LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Bflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorfluazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenoconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethamsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole; Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate E; Fenthion; Fenthion sulfone; Fenthion

sulfoxide; Fibronil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobil; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxyfop; Haloxyfop-2-Ethoxy-Ethyl(etotyl); Haloxyfop-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isafenphos; Isoproturon; Isoprazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaoxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyram; Norflurazon; Novaluron; Nuarimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Prometryne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Protham; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quazalofop-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Quazalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufonozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiadiazuron; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluand; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflururon; Triflurosulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide

285	Pestisit Analizi	KKMT/SOP01/PR15	Süt Ve Süt Ürünleri (Yem)	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
LC-MS/MS: 2,4-5-T; 2,4-D; 3,4,5-Trimethacarb; Acephate; Acetamiprid; Acetochlor; Acinbenzolar-S-Methyl; Alachlor; Aldicarb; Aldicarb Sulfone; Aldicarb Sulfoxide; Ametoctradin; Aminocarb; Amitraz; Asulam; Atrazine; Azamethiphos; Azimsulfuron; Azinphos-Ethyl; Azinphos-Methyl; Aziprotryne; Azoxystrobin; Beflubutamid; Benalaxyl; Bendiocarb; Benfuracarb; Bensulfuron-Methyl; Bentazone; Benthiavalicarb-Isopropyl; Benzoximate; Bifenazate; Bifenthrin; Bitertanol; Boscalid; Bromacil; Bromophos-Ethyl; Bromoxynil; Bromuconazole; Bupirimate; Buprofezin; Butoxycarboxim; Butralin; Butylate; Cadusafos; Carbaryl; Carbendazim+Benomyl; Carbofuran; Carbosulfan; Carboxin; Carfentrazone-Ethyl; Chlorantraniliprole; Chlorbromuron; Chlorfenvinphos; Chlorflazuron; Chloridazon; Chlorimuron-ethyl; Chloroxuron; Chlorpropham; Chlorpyrifos; Chlorpyrifos-Methyl; Chlorsulfuron; Chlorthiamid; Chromafenozide; Cinidon-ethyl; Clethodim; Clethodim sulfoxide; Clodinafop-propargyl+Clodinafop Propargyl Ester; Clofentezine; Cloquintocet-1-methylhexylester; Clothianidin; Cyanazine; Cyazofamid; Cycloate; Cyflufenamid; Cymoxanil; Cyproconazole; Cyprodinil; Cyromazine; Dazomet; Deltamethrin; Demeton-S-Methyl-Sulfone; Demeton-S-Methyl-Sulfoxide; Desmedipham; Diafenthiuron; Dialifos; Diallate; Diazinon; Dichlofluanid; Dichlorvos; Diclofop-Methyl; Dicrotophos; Diethofencarb; Difenconazole; Diflubenzuron; Dimefox; Dimethachlor; Dimethenamid; Dimethoate; Dimethomorph; Diniconazole; Dinitramine; Dinocap; Dioxacarb; Dioxathion; Diphenamid; Diphenylamine; Ditalimfos; Dithianon; Diuron; Amitraz met. DMF; Dodine; EPN; Epoxiconazole; EPTC; Ethametsulfuron-methyl; Ethiofencarb; Ethion; Ethiprole; Ethirimol; Ethofenprox; Ethofumasate; Ethoprophos; Ethoxyquin; Ethoxysulfuron; Etoazole; Famoxadone; Famphur; Fenamidone; Fenamiphos; Fenamiphos sulfone; Fenamiphos sulfoxide; Fenarimol; Fenazaquin; Fenbuconazole;						

Fenbutatin oxide; Fenhexamid; Fenoxaprop-p-ethyl; Fenoxycarb; Fenpropathrin; Fenpropidin; Fenpropimorph; Fenproximate; Fenpyroximate E; Fenthion; Fenthion sulfone; Fenthion sulfoxide; Fibronil; Fipronil sulfone; Flamprop-M-Isopropyl; Flazasulfuron; Florasulam; Fluazifop-P-butyl+Fluazifop-butyl; Fluazinam; Flubenzimine; Flucythrinate; Fludioxonil; Flufenoxuron; Flumioxazin; Fluometuron; Fluopicolide; Fluopyram; Fluoroglycofen-ethyl; Fluoxastrobilin; Fluquinconazole; Flurochloridone; Fluroxypyr; Flurtamone; Flusilazole; Flutolanil; Flutriafol; Fluxapyroxad; Fonofos; Forchlorfenuron; Formetanate; Fosthiazate; Fuberidazole; Furathiocarb; Halosulfuron methyl; Haloxypyr; Haloxypyr-2-Ethoxy-Ethyl(etotyl); Haloxypyr-R-methyl; Heptenophos; Hexaconazole; Hexaflumuron; Hexythiazox; Imazalil; Imazamox; Imazapic; Imazaquin; Imazethapyr; Imibenconazole; Imidacloprid; Indoxacarb; Iodosulfuron-Methyl-Sodium; Ioxynil; Ipconazole; Iprodione; Iprovalicarb; Isocarbophos; Isufenphos; Isoproturon; Isoprazam; Isoxaflutole; Kresoxim-Methyl; Lenacil; Linuron; Lufenuron; Malaaxon; Malathion; Mandipropamid; MCPA; Mecarbam; Mecoprop+Mecoprop-P; Mephosfolan; Mesosulfuron-methyl; Metalaxyl+Metalaxyl-M; Metamitron; Metconazole; Methacrifos; Methamidophos; Methamitron; Methidathion; Methiocarb; Methiocarb sulfone; Methiocarb sulfoxide; Methomyl; Methoxyfenozide; Metobromuron; Metolachlor+Metolachlor-S; Metolcarb; Metosulam; Metoxuron; Metrafenone; Metribuzin; Metsulfuron-methyl; Mevinphos; Molinate; Monocrotophos; Monolinuron; Monuron; Myclobutanil; Amitraz met. DMPF; Naled; Neburon; Nicosulfuron; Nitenpyr; Norflurazon; Novaluron; Nuaimol; Ofurace; Omethoate; Orthosulfamuron; Oxadiazon; Oxadixyl; Oxamyl; Oxasulfuron; Oxycarboxin; Paclobutrazol; Paraoxon-Ethyl; Paraoxon-methyl; Parathion-Ethyl; Pebulate; Penconazole; Pencycuron; Pendimethalin; Penoxsulam; Pethoxamid; Phenmedipham; Phenthoate; Phorate Sulfone; Phorate sulfoxide; Phosalone; Phosmet; Phosphamidon; Phoxim; Picolinafen; Piperonyl butoxide; Pirimicarb; Pirimicarb-desmethyl; Pirimiphos-ethyl; Pirimiphos-Methyl; Prochloraz; Profenofos; Profoxydim; Promecarb; Promeityne; Propamocarb; Propanil; Propaquizafop; Propargite; Propazine; Prophan; Propiconazole; Propoxur; Propyzamide; Proquinazid; Prosulfuron; Prothiofos; Pymetrozine; Pyraclostrobin; Pyrasulfotole; Pyrazophos; Pyridaben; Pyridaphenthion; Pyridate; Pyrimethanil; Pyriproxyfen; Quiazaflor-p-ethyl; Quinalphos; Quinclorac; Quinmerac; Quinomethionate (Chinomethionate); Quinoxifen; Quizalofop-p-ethyl; Rimsulfuron; Sethoxydim Solution; Simazine; Spinetoram (Spinosyn L+J); Spinosad (Spinosyn A+D); Spirodiclofen; Spirotetramat; Spiroxamine; Sulcotrione; Sulfotep; Tau-fluvalinate; Tebuconazole; Tebufenpyrad; Tebufenozide; Teflubenzuron; Tepraloxym; Terbutylazine; Terbutryn; Tetrachlorvinphos; Tetraconazole; Tetramethrin; Thiabendazole; Thiadiazuron; Thiamethoxam; Thidiazuron; Thifensulfuron Methyl; Thiobencarb; Thiodicarb; Thiofanox; Thiofanox-sulfoxide; Thiometon; Thiophanate-Methyl; Tolclofos-Methyl; Tolfenpyrad; Tolyfluanid; Topramezone; Tralkoxydim; Triadimefon; Triadimenol; Triallate; Triasulfuron; Triazophos; Tribenuron-methyl; Trichlorfon; Tricyclazole; Tridemorph; Trifloxystrobin; Triflumizole; Triflururon; Triflurosulfuron-methyl; Trinexapac-ethyl; Triticonazole; Tritosulfuron; Vamidothion; Zoxamide

286	Pestisit Analizi	KKMT/SOP01/PR15	Süt ve Süt Ürünleri	AOAC Official Method 2007.01 (Quechers)	Kromatografik	Kantitatif
GC-MS/MS: 1-Naphtylacetamide; 2,4-DDD; 2,4-DDE; 2,4-DDT; 2-Chloroaniline; 2-Phenylphenol; 3-Chloroaniline; 4,4-DDD; 4,4-DDE; 4,4-DDT; 4-Chloroaniline; Aclonifen; Acrinathrin, Benfluralin, Biphenyl, Bromocyclohexane, Bromopropylate, Captan, Carbophenothion, Chlorbenside, Chlorbufam, Chlordecone, Chlorfenapyr, Chlorfenson, Chlorobenzilate+Chloropropylate; Chloroneb, Chlorothalonil, Chlorthal-dimethyl, Chlozolate, Cyanophos, Cyfluthrin, Cyhalothrin-lambda, Cypermethrin, Demeton-S-Methyl, Dichlobenil, Dichlorobenzophenone-4,4, Dicloran, Dicofol, Dinobuton, Disulfoton, Endosulfan-alpha, Endosulfan-beta, Endosulfan-sulfate, Endrin, Endrin aldehyde, Endrin ketone, Esfenvalerate+Fenvalerate, Ethalfluralin, Etrimfos, Fenchlorphos, Fenitrothion, Fenson, Fluchloralin, Fluotrimazole, Flurprimidol, Folpet, Formothion, Halfenprox, HCH-alpha, HCH-beta, HCH-delta, HCH-gamma, Iodofenphos, Isocarbofos, Isodrin, Lactofen, Leptophos, Mefenpyr-diethyl, Methoxychlor, Mirex, Nitrpyrin, Nitrofen, Nitrothal-Isopropyl, Oxyfluorfen, Parathion-methyl, Pentachloroaniline, Permethrin, Perthane, Phorate, Procymidone, Profluralin, Quintozene, Resmethrin, Spiromesifen, Tecnazene, Tefluthrin, Terbacil, Terbufos, Tetradifon, Tetrasul, Trifluralin, Vinclozolin						

287	Pestisit Analizi Yüksek Polarlı	KKMT/SOP62/PR15	Gıda (Yüksek Yağ Ve Orta Düzeyde Su İçerikli Ürünler)	“Quick Method for the Analysis of Numerous Highly Polar Pesticides in Food Involving Extraction with Acidified Methanol and LCMS/MS Measurement I. Food of Plant Origin (QuPpe-PO-Method) Version 12”	Kromatografik	Kantitatif
LC-MS/MS: AMPA; Bromide Ion; Chlormequat Chloride; Cyanuric acid; Diquat dibromide monohydrate; Etephon-Hydroxy; Etephon; Fosetyl-Al; Glufosinate ammonium; Glyphosate; Maleic hydrazide; Mepiquat chloride; MPPA; N-Acetyl-AMPA; N-Acetyl-Glufosinate; Paraquat dichloride; Phosphonic acid; Sodium chlorate; Sodium perchlorate monohydrate						
288	Pestisit Analizi Yüksek Polarlı	KKMT/SOP62/PR15	Gıda (Yüksek Su İçerikli Ürünler; Yüksek Şeker ve Düşük Su İçerikli Ürünler; Yüksek Asit ve Yüksek Su İçerikli Ürünler; Yüksek Nişasta ve/veya Protein ve Düşük Su ve Yağ İçerikli Ürünler; Gıdalarda Yüksek Yağ ve Çok Düşük Su İçerikli Ürünler; Zor veya Özgün Ürünler)	“Quick Method for the Analysis of Numerous Highly Polar Pesticides in Food Involving Extraction with Acidified Methanol and LCMS/MS Measurement I. Food of Plant Origin (QuPpe-PO-Method) Version 12”	Kromatografik	Kantitatif
LC-MS/MS: AMPA; Bromide Ion; Chlormequat Chloride; Cyanuric acid; Diquat dibromide monohydrate; Etephon-Hydroxy; Etephon; Fosetyl-Al; Glufosinate ammonium; Glyphosate; Maleic hydrazide; Mepiquat chloride; MPPA; N-Acetyl-AMPA; N-Acetyl-Glufosinate; Paraquat dichloride; Phosphonic acid; Sodium chlorate; Sodium perchlorate monohydrate						
289	Pestisit Analizi Yüksek Polarlı	KKMT/SOP62/PR15	Hayvansal Kökenli Gıdalardan Elde Edilen Yağlar; Et (Kas) ve Su Ürünleri; Yumurta; Süt ve Süt Ürünleri	“Quick Method for the Analysis of Numerous Highly Polar Pesticides in Food Involving Extraction with Acidified Methanol and LCMS/MS Measurement II. Food of	Kromatografik	Kantitatif

				Animal Origin (QuPpe-AO-Method) Version 3.2”		
LC-MS/MS: AMPA; Bromide Ion; Chlormequat Chloride; Cyanuric acid; Diquat dibromide monohydrate; Etephon-Hydroxy; Etephon; Fosetyl-Al; Glufosinate ammonium; Glyphosate; Maleic hydrazide; Mepiquat chloride; MPPA; N-Acetyl-AMPA; N-Acetyl-Glufosinate; Paraquat dichloride; Phosphonic acid; Sodium chlorate; Sodium perchlorate monohydrate						
290	Pestisit Analizi Yüksek Polarlı	KKMT/SOP62/PR15	Yem (Yüksek Su İçerikli Ürünler; Yüksek Nişasta ve/veya Protein ve Düşük Su ve Yağ İçerikli Ürünler; Yüksek Yağ/Katı Yağ ve Düşük Su İçerikli Ürünler)	“Quick Method for the Analysis of Numerous Highly Polar Pesticides in Food Involving Extraction with Acidified Methanol and LCMS/MS Measurement I. Food of Plant Origin (QuPpe-PO-Method) Version 12”	Kromatografik	Kantitatif
LC-MS/MS: AMPA; Bromide Ion; Chlormequat Chloride; Cyanuric acid; Diquat dibromide monohydrate; Etephon-Hydroxy; Etephon; Fosetyl-Al; Glufosinate ammonium; Glyphosate; Maleic hydrazide; Mepiquat chloride; MPPA; N-Acetyl-AMPA; N-Acetyl-Glufosinate; Paraquat dichloride; Phosphonic acid; Sodium chlorate; Sodium perchlorate monohydrate						
291	Fosfin Analizi	KKMT/SOP63/PR15	Gıdalar	EURL Analysis of the Fumigant Phosphine	GC-MS MS Headspace	Kantitatif
292	Osmofilik Maya Sayımı	MAL/SOP65/PR15	Meyve ve Şekerli Ürünler	ISO 21527-2	Yayma Plak	Kantitatif
293	<i>Escherichia coli</i> EMS Yöntemi ile Tespiti ve Sayımı	MAL/SOP66/PR15	Gıda (Süt ve Süt Ürünleri Hariç), Yem Ve Gıda İşletmesinden Alınan Çevresel Numuneler	ISO 16649-3 ISO 18593	EMS yöntemi	Kantitatif/ Kalitatif
294	GHB811 Pamuk Tip Belirleme Analizi ^{1 2}	MBAL/SOP95/PR15	Gıda ve Yem	EURL Metot (QT-EVE-GH-013), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
295	DP-073496-4 Kanola Tip Belirleme Analizi ^{1 2}	MBAL/SOP96/PR15	Gıda ve Yem	EURL Metot (QT-EVE-BN-009), ISO 21569, ISO 24276	Real Time PCR	Kalitatif

296	GDO Tarama Analizi (p35S, tNOS, pFMV) ^{1 2}	MBAL/SOP01/PR15	Gıda üretimi amacıyla yetiştirilen ürünlerin yaprak, gövde, kök gibi bitki kısımları	Kit Metodu (Eurofins GeneScan GMO Screen RT 35S/NOS/FMV IPC Test Kit Manual), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
297	GDO Tarama Analizi (p35S, tNOS,pFMV) ^{1 2}	MBAL/SOP56/PR15	Gıda üretimi amacıyla yetiştirilen ürünlerin yaprak, gövde, kök gibi bitki kısımları	Kit Metodu (Bosphore® GMO Screening Kit S35-TNOS- FMV), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
298	GDO Tarama Analizi (p35S, tNOS) ^{1 2}	MBAL/SOP82/PR15	Gıda üretimi amacıyla yetiştirilen ürünlerin yaprak, gövde, kök gibi bitki kısımları	EURL Metot (QL-ELE-00-013), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
299	GDO Tarama Analizi (pFMV) ^{1 2}	MBAL/SOP83/PR15	Gıda üretimi amacıyla yetiştirilen ürünlerin yaprak, gövde, kök gibi bitki kısımları	EURL Metot (QL-ELE-00-015), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
300	Ani Sıcaklık Değişimine Dayanım	FAL/SOP51/PR15	Cam Esaslı Ambalaj	TS 2913	Organoleptik	Kalitatif
301	Boya Geçirgenliği Analizi (sulu, asitli ve yağlı)	FAL/SOP52/PR15	Gıda İle Temas Eden Baskılı Kağıt Karton Ve Plastik Madde Ve Malzemeler	EN 646, EN 645	Organoleptik	Kalitatif
302	Yabancı Madde Muhtevası Tayini	FAL/SOP11/PR15	Yağlı Tohumlar	TS 2947 EN ISO 658	Gravimetrik	Kantitatif
303	Ham Protein Tayini – Dumas	KAL/SOP108/PR15	Gıda ve Yem	ISO 16634-1 ve ISO 16634-2	Dumas Metot	Kantitatif
304	Ham Protein Tayini - Kjeldahl	KAL/SOP109/PR15	Yem	ISO 5983-2	Titrimetrik Kjeldahl	Kantitatif
305	Ham Protein Tayini - Kjeldahl	KAL/SOP110/PR15	Hububatlar	ISO 20483	Titrimetrik Kjeldahl	Kantitatif
306	Tersiyer Bütıl Alkol (TBA) Analizi	KAL/SOP111/PR15	Alkollü İçkiler	İşletme İçi Metot- 'KAL/SOP111/PR15' (CONSLEG:2000R2870 ve TKG 2017/9 dan modifiye edilmiştir.)	Kromatografik GC-FID	Kantitatif

307	Hidroksi Metil Furfural (HMF) Tayini	KKMT/SOP69/PR15	Meyve-Sebze Mamulleri, Bal, Reçel, Pekmez	TS 6178 ISO 7466, IHC Bölüm 5.1	Kromatografik HPLC-UV/DAD	Kantitatif
308	Laktoz Tayini	KKMT/SOP70/PR15	Süt ve Süt Ürünleri	ISO 22662	Kromatografik HPLC-RID	Kantitatif
309	Şeker Bileşenleri Tayini (Glikoz, Fruktoz, Sakkaroz, Maltoz, Laktoz)	KKMT/SOP71/PR15	Meyve-Sebze Mamulleri, Bal, Reçel, Pekmez	TS 13359, TS EN 12630 IHC Bölüm 7.2	Kromatografik HPLC-RID	Kantitatif
310	Benzoil Peroksit Tayini	KKMT/SOP72/PR15	Un	Shimadzu Aplikasyon Notu	Kromatografik HPLC-UV/DAD	Kantitatif
311	Propineb Tayini	KKMT/SOP73/PR15	Bebek formülleri, devam formülleri ve bebek ve küçük çocuk ek gıdaları	İşletme İçi Metot- 'KKMT/SOP73/PR15' (A Rapid and Sensitive Analysis of Dithiocarbamate Fungicides Using Modified QuEChERS Method and Liquid Chromatography-Tandem Mass Spectrometry'den modifiye edilmiştir.)	Kromatografik LC-MS/MS	Kantitatif
312	Veteriner İlaç Kalıntıları Tayini	KKMT/SOP19/PR15	Kırmızı Et Kas Doku (Sığır, Koyun ve Keçi)	Journal of Chromatography A, 1213(2008) 189-199 Analytica Chimica Acta 637 (2009) 87-91 Journal of Chromatography A, 1216(2009) 8217-8223 Analytica Chimica Acta 483 (2003) 99-109 Analytica Chimica Acta 880 (2015) 103-121 Journal of Chromatography A, 1217 (2010) 6394-6404	Kromatografik LC-MS/MS	Kantitatif

LC-MS/MS (53 Etken Madde)

Antihelmintikler Grubu-3 Adet

Levamisol; Thiabendazole 5 Hydroxy; Thiabendazole

Quinolone & Fluoroquinolone Grubu -12 Adet

Nalidixic Acid; Sarafloxacin; Enoxacin ; Lomefloxacin; Ofloxacin ; Siprofloksasin; Danofloksasin ; Difloksasin ; Enrofloksasin ; Marbofloksasin ; Flumequin ; Oksolinik asit

Macrolit Grubu-5 Adet

Spiramycin; Neospiramycin; Eritromisin ; Tilmikosin ; Tilosin ;

Sulfanamid Grubu-12 Adet

Sulfabenzamide; Sulfadoxine; Sulfameter; Sulfisoxazole ; Sülfakloropiridazin, Sülfadiazine, Sülfadimetoksin, Sülfamerazine, Sülfametazin, Sülfametaksazol, Sülfakuinoksalin, Sülfatiazol

Tetrasiklin Grubu-7 Adet

Chlortetracycline; Oxytetracycline; Tetracycline; 4-epi-Tetracycline; 4-epi-Oxytetracycline; 4-epi-Chlortetracycline ; Doksisiklin

Penisilin Grubu – 7 Adet

Amoksisilin, Ampisilin, Kloksasilin, Dikloksasilin, Nafsilin, Oksasilin, Pen-G(Benzilpenisilin)

Diaminoprimidin Grubu-1 Adet

Trimethoprim

Fenikol Grubu – 4 Adet

Florfenikol, Tiamfenikol, Kloramfenikol, Florfenikol amine

Beta-Lactam Grubu – 2 Adet

Sefkuinom ; Sefaleksim,

313	Veteriner İlaç Kalıntıları Tayini	KKMT/SOP19/PR15	Kanatlı Kas Doku (Tavuk, Hindi)	Journal of Chromatography A, 1213(2008) 189-199 Analytica Chimica Acta 637 (2009) 87-91 Journal of Chromatography A, 1216(2009) 8217-8223 Analytica Chimica Acta	Kromatografik LC-MS/MS	Kantitatif
-----	--------------------------------------	-----------------	------------------------------------	--	---------------------------	------------

				483 (2003) 99–109 Analytica Chimica Acta 880 (2015) 103–121 Journal of Chromatography A, 1217 (2010) 6394–6404		
LC-MS/MS (53 Etken Madde) <u>Antihelmintikler Grubu-3 Adet</u> Levamisol; Thiabendazole 5 Hydroxy; Thiabendazole <u>Quinolone & Fluoroquinolone Grubu-13 Adet</u> Nalidixic Acid; Sarafloxacin; Enoxacin; Lomefloxacin; Ofloxacin; Marbofloksasin ; Siprofloksasin; Danofloksasin ; Difloksasin ; Enrofloksasin ; Marbofloksasin ; Flumequin ; Oksolinik asit <u>Macrolit Grubu- 5 Adet</u> Spiramycin; Neospiramycin ; Eritromisin ; Tilmikosin ; Tilosin <u>Sulfanamid Grubu-12 Adet</u> Sulfabenzamide; Sulfadoxine; Sulfameter; Sulfisoxazole ; Sülfakloropiridazin, Sülfadiazine, Sülfadimetoksin, Sülfamerazine, Sülfametazin, Sülfametaksozol, Sülfakuinoksalin, Sülfatiazol <u>Tetrasiklin Grubu-7 Adet</u> Chlortetracycline; Oxytetracycline; Tetracycline; 4-epi-Tetracycline; 4-epi-Oxytetracycline; 4-epi-Chlortetracycline; Doksisiklin ; <u>Penisilin Grubu -7 adet</u> Amoksisilin, Ampisilin, Kloksasilin, Dikloksasilin, Nafsilin, Oksasilin, Pen-G(Benzilpenisilin) <u>Beta-Laktam Grubu -1 adet</u> Cefalexin <u>Diaminoprimidin Grubu-1 Adet</u> Trimethoprim <u>Fenikol Grubu - 4 Adet</u> Florfenikol, Tiamfenikol, Kloramfenikol, Florfenikol amine						

314	Veteriner İlaç Kalıntıları Tayini	KKMT/SOP19/PR15	Balık Kas Doku	Journal of Chromatography A, 1213(2008) 189-199 Analytica Chimica Acta 637 (2009) 87-91 Journal of Chromatography A, 1216(2009) 8217-8223 Analytica Chimica Acta 483 (2003) 99-109 Analytica Chimica Acta 880 (2015) 103-121 Journal of Chromatography A, 1217 (2010) 6394-6404	Kromatografik LC-MS/MS	Kantitatif
LC-MS/MS (53 Etken Madde) <u>Antihelmintikler Grubu-3 Adet</u> Levamisol; Thiabendazole 5 Hydroxy; Thiabendazole <u>Quinolone & Fluoroquinolone Grubu -13 Adet</u> Nalidixic Acid; Sarafloxacin; Enoxacin; Lomefloxacin; Ofloxacin; Marbofloksasin ; Siprofloksasin; Danofloksasin ; Difloksasin ; Enrofloksasin ; Marbofloksasin ; Flumequin ; Oksolinik asit <u>Macrolit Grubu-5 Adet</u> Spiramycin; Neospiramycin ; Eritromisin ; Tilmikosin ; Tilosin <u>Sulfanamid Grubu-12 Adet</u> Sulfabenzamide; Sulfadoxine; Sulfameter; Sulfisoxazole ; Sülfakloropiridazin, Sülfadiazine, Sülfadimetoksin, Sülfamerazine, Sülfametazin, Sülfametaksozol, Sülfakuinoksalin, Sülfatiazol <u>Tetrasiklin Grubu-7 Adet</u> Chlortetracycline; Oxytetracycline; Tetracycline; 4-epi-Tetracycline; 4-epi-Oxytetracycline; 4-epi-Chlortetracycline; Doksisiklin <u>Penisilin Grubu -7 adet</u> Amoksisilin, Ampisilin, Kloksasilin, Dikloksasilin, Nafsilin, Oksasilin, Pen-G(Benzilpenisilin) <u>Beta-Laktam Grubu -1 adet</u> Cefalexin						

<u>Diaminoprimidin Grubu-1 Adet</u> Trimethoprim <u>Fenikol Grubu - 4 Adet</u> Florfenikol, Tiamfenikol, Kloramfenikol, Florfenikol amine						
315	Veteriner İlaç Kalıntıları Tayini	KKMT/SOP19/PR15	Süt	Journal of Chromatography A, 1213(2008) 189-199 Analytica Chimica Acta 637 (2009) 87-91 Journal of Chromatography A, 1216(2009) 8217-8223 Analytica Chimica Acta 483 (2003) 99-109 Analytica Chimica Acta 880 (2015) 103-121 Journal of Chromatography A, 1217 (2010) 6394-6404	Kromatografik LC-MS/MS	Kantitatif
LC-MS/MS (48 Etken Madde) <u>Antihelmintikler Grubu-3 Adet</u> Levamisol; Thiabendazole 5 Hydroxy; Thiabendazole <u>Quinolene & Fluoroquinolone Grubu -9 Adet</u> Nalidixic Acid; Sarafloxacin; Enoxacin; Lomefloxacin; Ofloxacin; Marbofloksasin; Ciprofloxacin; Danofloksasin ; Flumequin <u>Macrolit Grubu-5 Adet</u>						

Spiramycin; Neospiramycin ; Eritromisin ; Tilmikosin ; Tilosin

Sulfanamid Grubu-12 Adet

Sulfabenzamide; Sulfadoxine; Sulfameter; Sulfisoxazole ; Sülfakloropiridazin, Sülfadiazine, Sülfadimetoksin, Sülfamerazine, Sülfametazin, Sülfametaksazol, Sülfakuinoksalin, Sülfatiazol

Penisilin Grubu -7 adet

Amoksisilin, Ampisilin, Kloksasilin, Dikloksasilin, Nafsilin, Oksasilin, Pen-G(Benzilpenisilin)

Tetrasiklin Grubu-7 Adet

Chlortetracycline; Oxytetracycline; Tetracycline; 4-epi-Tetracycline; 4-epi-Oxytetracycline; 4-epi-Chlortetracycline; Doksisiklin

Diaminoprimidin Grubu-1 Adet

Trimethoprim

Beta-Lactam Grubu – 2 Adet

Sefkuinom ; Sefaleksine,

Fenikol Grubu – 2 Adet

Tiamfenikol, Kloramfenikol

316	Veteriner İlaç Kalıntıları	KKMT/SOP19/PR15	Tereyağı	Journal of Chromatography A, 1213(2008) 189-199 Analytica Chimica Acta 637 (2009) 87-91 Journal of Chromatography A, 1216(2009) 8217-8223 Analytica Chimica Acta 483 (2003) 99–109 Analytica Chimica Acta 880 (2015) 103–121 Journal of Chromatography A, 1217 (2010) 6394–6404	Kromatografik LC-MS/MS	Kantitatif
LC-MS/MS (53 Etken Madde) <u>Antihelmintikler Grubu-3 Adet :</u> Levamisol; Thiabendazole 5 Hydroxy; Thiabendazole						

Quinolone & Fluoroquinolone Grubu -12 Adet

Enoxacin; Lomefloxacin; Ofloxacin ; Siprofloksasin, Danofloksasin, Difloksasin, Enrofloksasin, Flumequin, Marbofloksasin, Nalidiksik asit, Oksolinik asit, Sarafloksasin

Macrolit Grubu-5 Adet

Spiramycin; Neospiramycin ; Eritromisin ; Tilmikosin ; Tilosin

Sulfanamid Grubu-12 Adet

Sulfabenzamide; Sulfadoxine; Sulfameter; Sulfisoxazole ; Sülfakloropiridazin, Sülfadiazine, Sülfadimetoksin, Sülfamerazine, Sülfametazin, Sülfametaksazol, Sülfakuinoksalin, Sülfatiazol

Tetrasiklin Grubu-7 Adet

Chlortetracycline; Oxytetracycline; Tetracycline; 4-epi-Tetracycline; 4-epi-Oxytetracycline; 4-epi-Chlortetracycline ; Doksisisiklin

Penisilin Grubu -8 adet

Cephapirin, Amoksisilin, Ampisilin, Kloksasilin, Dikloksasilin, Nafsilin, Oksasilin, Pen-G(Benzilpenisilin)

Beta-Laktam Grubu -3 adet

Ceftiofur HCl ; Sefkuinom ; Sefaleksim,

Diaminoprimidin Grubu-1 Adet

Trimethoprim

Fenikol Grubu – 4 Adet

Florfenikol, Tiamfenikol, Kloramfenikol, Florfenikol amine

317	Veteriner İlaç Kalıntıları	KKMT/SOP19/PR15	Yem	Journal of Chromatography A, 1213(2008) 189-199 Analytica Chimica Acta 637 (2009) 87-91 Journal of Chromatography A, 1216(2009) 8217-8223 Analytica Chimica Acta 483 (2003) 99–109 Analytica Chimica Acta 880 (2015) 103–121	Kromatografik LC-MS/MS	Kantitatif
-----	----------------------------	-----------------	-----	--	------------------------	------------

				Journal of Chromatography A, 1217 (2010) 6394–6404		
LC-MS/MS (53 Etken Madde) Antihelmintikler Grubu-3 Adet: Levamisol; Thiabendazole 5 Hydroxy; Thiabendazole Quinolene & Fluoroquinolone Grubu -12 Adet Enoxacin ; Lomefloxacin; Ofloxacin ; Sipprofloksasin, Danofloksasin, Difloksasin, Enrofloksasin, Flumekuın, Marbofloksasin, Nalidiksik asit, Oksolinik asit, Sarafloksasin Macrolit Grubu-5 Adet Spiramycin; Neospiramycin ; Eritromisin ; Tilmikosin ; Tilosin Nitroimidazol Grubu-1 Adet Ronidazole Sulfanamid Grubu-12 Adet Sulfabenzamide; Sulfadoxine; Sulfameter; Sulfisoxazole ; Sülfakloropiridazin, Sülfadiazine, Sülfadimetoksin, Sülfamerazine, Sülfametazin, Sülfametaksozol, Sülfakuinoksalin, Sülfatiazol, Tetrasiklin Grubu-7 Adet Chlortetracycline; Oxytetracycline; Tetracycline; 4-epi-Tetracycline; 4-epi-Oxytetracycline; 4-epi-Chlortetracycline; Doksisisiklin Penisilin Grubu -8 adet Cephapirin ; Amoksisilin, Ampisilin, Kloksasilin, Dikloksasilin, Nafsilin, Oksasilin, Pen-G(Benzilpenisilin) Beta-Laktam Grubu -3 adet Ceftiofur HCl ; Sefkuinom ; Sefaleksın, Fenikol Grubu – 4 Adet Florfenikol, Tiamfenikol, Kloramfenikol, Florfenikol amine Diaminoprimidin Grubu-1 Adet Trimethoprim						
318	MON 94100 Kanola Tip Belirleme Analizi ^{1,2}	MBAL/SOP98/PR15	Gıda ve Yem	EURL Metot (QT-EVE- BN-014), ISO 21569, ISO 24276	Real Time PCR	Kalitatif
319	Sülfat İndirgeyen Clostridium türlerinin Koloni Sayımı	MAL/SOP67/PR15	Gıda	ISO 15213-1 ISO 18593	Dökme Plak	Kantitatif

320	Kurşun (Pb) Analizi	KKMT/SOP76/PR15	Şarap	OIV-OENO 666-2023	Spektrofotometrik ICP-MS	Kantitatif
321	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A)	KKMT/SOP03/PR15	Yağlı Tohumlar ve ürünleri	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
322	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A)	KKMT/SOP03/PR15	Sert Kabuklu Meyveler ve ürünleri	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
323	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A)	KKMT/SOP03/PR15	Kurutulmuş Meyveler ve ürünleri	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
324	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A)	KKMT/SOP03/PR15	Kurutulmuş Sebzeler ve ürünleri	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
325	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A, Zearalenon, Deoksinivalenol, Fumonisin B1, Fumonisin B2, Fumonisin B1+B2, T2-HT2)	KKMT/SOP03/PR15	Tahıl ve Tahıl Ürünleri	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif

326	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A)	KKMT/SOP03/PR15	Baharatlar	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
327	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A, Zearalenon, Deoksinivalenol, Fumonisin B1, Fumonisin B2, Fumonisin B1+B2)	KKMT/SOP03/PR15	Bebek ve küçük çocuk Ek Gıdaları	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
328	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A, Zearalenon, Deoksinivalenol, Fumonisin B1, Fumonisin B2, Fumonisin B1+B2, T2- HT2)	KKMT/SOP03/PR15	Yem, Yem Katkı ve Yem Maddeleri	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
329	Aflatoksin M1	KKMT/SOP04/PR15	Süt ve Süt Tozu	İşletme içi metot KKMT/SOP04/PR15 (Aflatoksin M1 Jasem Method)	Kromatografik	Kantitatif
330	Aflatoksin M1	KKMT/SOP04/PR15	Bebek ve Devam Formülleri/Bebek sütü ve devam sütü	İşletme içi metot KKMT/SOP04/PR15 (Aflatoksin M1 Jasem Method)	Kromatografik	Kantitatif

331	Aflatoksin M1	KKMT/SOP04/PR15	Süt ürünleri	İşletme içi metot KKMT/SOP04/PR15 (Aflatoksin M1 Jasem Method)	Kromatografik	Kantitatif
332	Okratoksin-A	KKMT/SOP03/PR15	Şarap, aromatize şarap bazlı içki ve aromatize şarap kokteyli, Bira, malt içeceği, Alkollü içecekler üzüm suyu	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
333	Okratoksin-A	KKMT/SOP03/PR15	Çiğ Kahve	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
334	Okratoksin-A	KKMT/SOP03/PR15	Kavrulmuş kahve, suda çözünebilen kahve ve ekstraktları	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
335	Okratoksin-A	KKMT/SOP03/PR15	Kakao ve Kakao ürünlerinde	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
336	Okratoksin-A	KKMT/SOP03/PR15	Kurutulmuş bitkiler (baharat ve kuru meyveler hariç)	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
337	Patulin	KKMT/SOP05/PR15	Meyve suları, konsantreden elde edilen meyve suları ve meyve nektarları, meyveli içecekler (TGK Meyve Suyu ve Benzeri Ürünler Tebliği kapsamındaki ürünler)	TS EN 14177	Kromatografik	Kantitatif
338	Patulin	KKMT/SOP77/PR15	Katı elma ürünleri	İşletme içi metot KKMT/SOP77/PR15 (TS EN 14177' den modifiye edilmiştir.)	Kromatografik	Kantitatif
339	Patulin	KKMT/SOP77/PR15	Distile Alkollü İçecekler, elma şarabı ve elmadan elde edilen fermente içecekler, diğer alkollü içkiler	İşletme içi metot KKMT/SOP77/PR15 (TS EN 14177' den modifiye edilmiştir.)	Kromatografik	Kantitatif

340	Patulin	KKMT/SOP77/PR15	Bebek ve küçük çocuklar için elma ürünleri ve ek gıdalar	İşletme içi metot KKMT/SOP077/PR15 (TS EN 14177' den modifiye edilmiştir.)	Kromatografik	Kantitatif
341	Fumonisin B1,B2	KKMT/SOP03/PR15	Mısır ve mısırdan üretilen gıda ürünleri	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
342	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A, Zearalenon, Deoksinivalenol, Fumonisin B1, Fumonisin B2, Fumonisin B1+B2, T2-HT2)	KKMT/SOP03/PR15	Bitkisel Çaylar (Siyah çay ve yeşil çay hariç)	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
343	Multitoksin (Aflatoksin B1-B2-G1-G2, Toplam Aflatoksin B1+B2+G1+G2, Okratoksin A, Zearalenon, Deoksinivalenol, Fumonisin B1, Fumonisin B2, Fumonisin B1+B2, T2-HT2)	KKMT/SOP03/PR15	Bitkisel Sıvı Yağlar	İşletme içi metot KKMT/SOP03/PR15 (Multitoksin Jasem Method)	Kromatografik	Kantitatif
344	Dioksin, Furan, dPCBs ve İndikatör PCB Analizi	KKMT/SOP64/PR15	Gıda, Gıda Katkıları, Gıda Takviyeleri, Bebek ve Küçük Çocuk Ek Gıdalarında (Yaş Ağırlık Üzerinden Verilen)	EPA1613B EPA 1668C	Kromatografik (GC-MS/MS)	Kantitatif

345	Dioksin, Furan, dPCBs ve İndikatör PCB Analizi	KKMT/SOP65/PR15	Yem, Yem Hammaddeleri, Yem Katkı Maddeleri ve Yemlik Yağ ve Ürünlerinde	TS EN 16215	Kromatografik (GC-MS/MS)	Kantitatif
346	Dioksin, Furan, dPCBs ve İndikatör PCB Analizi	KKMT/SOP66/PR15	Gıda, Gıda Katkıları, Gıda Takviyelerinde (Gram Yağ Üzerinden Verilen)	EPA1613B EPA 1668C	Kromatografik (GC-MS/MS)	Kantitatif
347	Polisiklik Aromatik Hidrokarbonların (PAH) [Benzo(a)pyren, Krisen, Benzo(a)antrasen, Benzo(b)floranten] Tayini	KKMT/SOP67/PR15	Tüm Gıdalar, Bebek ve Küçük Çocuk Gıdası, Bebek Formülleri ve Devam Formülleri	İşletme İçi Metot: “KKMT/SOP67/PR15” (Aplikasyon notu: 30126 yayınından modifiye edilmiştir) (AOAC Official Method 2007.01 standardından modifiye edilmiştir)	Kromatografik (GC-MS/MS)	Kantitatif

*Kullanılan Teknik/Ölçüm Tekniği: Spektrofotometrik, Kromatografik (GC, GC-MS, GC-MS/MS, LC-MS, LC-MS/MS, HPLC v.b), Gravimetrik, Volumetrik, Elektroforez, Organoleptik, Mikroskopik vb.

¹ Bu analizlere ilişkin çalışma izni gıda ve yem ürünlerinde **DNA ekstraksiyonunu** [Eurofins GENESpin DNA İzolasyon Kit Metodu, ISO 21571] içermektedir.

² Bu analizlere ilişkin çalışma izni gıda, yem, jelatin ve jelatin içerikli ürünlerinde **DNA ekstraksiyonunu** [EurofinsGeneScan-GENESpin ve QIAGEN DNeasy Mericon Food, ISO 21571] içermektedir.